



 *design and quality*
MADE IN ITALY since 1980

O R M E S A®



M E A S U R E S



BIRILLO 1




BIRILLO 2



BIRILLO 3



BIRILLO 4



THE NUMBERS DIVIDED BY THE DASH SPECIFY A MINIMUM AND A MAXIMUM ADJUSTMENT

Sizes	Max. user mass	Weight of the aid	Turning width	A	B	C	D	E	F	G
1	kg 30	kg 8,5	cm 67	cm 18	cm 59	cm 50	cm 19-27	cm 36-50	cm 18	cm 22
2	kg 40	kg 10	cm 73	cm 19	cm 65	cm 51	cm 22-30	cm 40-56	cm 20	cm 26
3	kg 50	kg 11	cm 81	cm 21	cm 73	cm 59	cm 22-34	cm 42-62	cm 23	cm 30
4	kg 60	kg 12,8	cm 89	cm 28	cm 78	cm 60	cm 22-40	cm 50-76	cm 30	cm 38

A C C E S S O R I E S



860 - concave headrest
adjustable in height



824 - transparent table

BIRILLO safeness is in compliance with the **93/42 European Directive**:

- The **frame**, adjustable in height, is in steel and it is painted with epoxy-atossic varnish;
- Fire-retardant **upholstery**;
- The **seat** is padded and covered with plastic material, resisting to wear and tear and cleaning;
- The elastic **wall-protector buffers** are over the wheels;
- The **solid tyres** are castored and they run on ball bearings;
- The **handlebar** is adjustable in height and depth.

BIRILLO was conceived and thought to let children from 2 to 8 years old precociously move by themselves in the sitting position. **Its upholstery's fire-retardant. The triangular sitting position** prevents from the extreme hip and trunk extension and lets the child have his head and trunk well lined up. According to this, the aid helps the child having a **right posture** and his **moving impulse**. **The abduction block** keeps the pelvis in the triangular sitting. **Its posture stimulates the movement** of the arms and the rotation of the trunk. Thanks to its **shape** and **colours**, **BIRILLO** is **accepted as a "toy"** by the child and this predisposes to a serene use, without fear, by **stimulating his curiosity of exploring the surrounding environment**. For this reason **BIRILLO** plays a very important part in the **development of the perceptive and motor visual functions of the child**, which are necessary to know and represent the space.

TESTED BY **TÜV** BERLIN

FDA REGISTERED (USA)

CE MARKED

Authorized Dealer

ORMESA® s.r.l.
Via A. Da Sangallo, 1 - PACIANA
06034 FOLIGNO (PG) ITALY
TEL. +39 0742 22927 FAX +39 0742 22637
E-mail: info@ormesa.com - www.ormesa.com

COMPANY WITH QUALITY MANAGEMENT
SYSTEM CERTIFIED BY DNV
= **ISO 9001:2008** =