

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



Certified by **CERAH France**
€ marked


- F**UNCTIONS
- Steel and aluminum folding **frame** with a non-toxic epoxy varnish finish
 - Fire-retardant padding and removable **upholstery**
 - Rigid **backrest** with 4 quick release recline positions
 - Fixed rigid **seat**
 - Adjustable and tip-up **footrest**
 - Drum **brakes**
 - **Wheels:** front puncture proof cm 17,5 castor wheels, and back pneumatic cm 30 wheels, acting as dampers.

Please check off the interested products:

MEASURES

- OBI-2**
- OBI-3**
- OBI-4**

Frame colour: **blue**
Upholstery colour: **green-pink-blue** (combined)



SIZE	max load cap	weight of the aid	A	B	C	D	E	G	H	I	L	M
2	Kg. 50	Kg. 17,4	cm 62	cm 33	cm 30	cm 23-35	18°	cm 60	cm 97	cm 99	cm 92	cm 50
3	Kg. 60	Kg. 18,5	cm 68	cm 37	cm 36	cm 29-44	18°	cm 62,5	cm 103	cm 109	cm 95	cm 46
4	Kg. 70	Kg. 20,1	cm 76	cm 41	cm 40	cm 32-47	18°	cm 66,5	cm 105	cm 124	cm 103	cm 51

ACCESORIES

<input type="checkbox"/> 812 Wheel direction locks 	<input type="checkbox"/> 816 Side supports 	<input type="checkbox"/> 817 Leg straps 	<input type="checkbox"/> 818 Thermal cover 	<input type="checkbox"/> 819 Canopy 	<input type="checkbox"/> 824 Table with edge and hollow 	
<input type="checkbox"/> 825 Rain cover (it can only be attached after the canopy 819) 	<input type="checkbox"/> 827 Foot straps 	<input type="checkbox"/> 828 Seat belt 	<input type="checkbox"/> 832 Optional frame with four wheels with solid tyres; cm 20 rear wheels, cm 17 front wheels 	<input type="checkbox"/> 833 Optional frame with four cm 30 wheels 	<input type="checkbox"/> 834 Abduction block 	<input type="checkbox"/> 835 Headrest with two side elements
<input type="checkbox"/> 837 Tray 	<input type="checkbox"/> 839 Front handle 	<input type="checkbox"/> 853 Adjustable harness 	<input type="checkbox"/> 858 Basket 	<input type="checkbox"/> 891 Set of 4 tie down hooks for the forward-facing transport of the pushchair in a motor vehicle (private vans/cars, buses...) 	<input type="checkbox"/> 894 45° Pelvic belt 	

Notes:

