




MARKED 

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

FEATURES

- **Frame** height-adjustable, made of steel painted with non-toxic epoxy paint.
- Fire-retardant **upholstery**.
- **Seat** triangular to hold the pelvis (seat/back at a 90° angle). It is padded and wrap-around, upholstered with strong, washable plastic.
- **Bumpers** elastic wall-protector bumpers and protective balls on the upper ends of the rods.
- **Wheels** 80 mm fully-swivelling, rolling on ball-bearings, 2 with brakes.
- **Handle bar** adjustable in both height and depth.
- **Abduction block** adjustable to contain the pelvis.



Please check off the interested products



MEASURES

BIRILLO 1 Extra Small





Frame colour: **blue**

Upholstery colours: **yellow-green-red**

BIRILLO 2 Small

BIRILLO 3 Medium

BIRILLO 4 Large





SIZE	Aid Weight	Turning Width	Max. Cap. Load	A	B	C	D	E	F	G	H
1	Kg. 8,5	cm 67	Kg. 30	cm 18	cm 59	cm 50	cm 20-28	cm 34-48	cm 18	cm 22	cm 70-95
2	Kg. 10	cm 73	Kg. 40	cm 19	cm 65	cm 51	cm 23-31	cm 40-56	cm 20	cm 26	cm 85-105
3	Kg. 11	cm 81	Kg. 50	cm 21	cm 73	cm 59	cm 23-35	cm 40-65	cm 23	cm 30	cm 95-120
4	Kg. 12,8	cm 89	Kg. 60	cm 26	cm 78	cm 60	cm 22-40	cm 43-74	cm 30	cm 38	cm 110-140

ADDITIONAL COMPONENTS

812 Directional lock

824 Transparent tray

860 Headrest



Notes:

