

- Tested by **BERLIN CERT GmbH** and listed in the **GERMAN HMV** (Technical Aids Listing)
- Certified by **CERAH France**
- **FDA** registered (USA)
- **CE** Marked

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

**mis. Small**

Upholstery colour:  **dark grey/green**

**mis. Medium**

**dark grey/pink**

Frame colour: **black**

**light grey/green**

**light grey/pink**



**Attention!**  
Read carefully the user and maintenance manual

Hypen-separated numbers indicate the minimum and maximum range of adjustment

MEASURES



Bug with 856 indoor base

	A	B	C	D	E	F	G	J	H	I	L	M	N	O	P	max. load	frame weight	seat weight	total weight
<b>Small</b>	20-30 cm 7,8-11,8 in	52-62 cm 20,4-24,4 in	37 cm 14,5 in	47-75 cm 18,5-29,5 in	20-30 cm 7,8-11,8 in	14-30 cm 5,5-11,8 in	0°-35°	35°	95-132 cm 37,4-51,9 in	56 cm 22 in	34-61 cm 13,3-24 in	80 cm 31,4 in	40 cm 15,7 in	25 cm 9,8 in	90-120 cm 35,4-47,2 in	35 kg 77,1 lb	14,5 kg 31,9 lb	8,5 kg 18,7 lb	23 kg 50,7 lb
<b>Medium</b>	25-36 cm 9,8-14,1 in	60-75 cm 23,6-29,5 in	46 cm 18,1 in	47-75 cm 18,5-29,5 in	30-40 cm 11,8-15,7 in	20-38 cm 7,8-14,9 in	0°-25°	40°	106-148 cm 41,7-58,2 in	62 cm 24,4 in	34-61 cm 13,3-24 in	83 cm 32,6 in	45 cm 17,7 in	32 cm 12,5 in	115-140 cm 45,2-55,1 in	45 kg 99,2 lb	16,2 kg 35,7 lb	12,7 kg 27,9 lb	28,9 kg 63,7 lb



Bug with 851 indoor base

	A	B	C	D	E	F	G	J	H	I	L	M	N	O	P	max. load	frame weight	seat weight	total weight
<b>Small</b>	20-30 cm 7,8-11,8 in	52-62 cm 20,4-24,4 in	37 cm 14,5 in	34-50 cm 13,3-19,6 in	20-30 cm 7,8-11,8 in	14-30 cm 5,5-11,8 in	0°-35°	35°	84-110 cm 33-43,3 in	50 cm 19,6 in	23-39 cm 9-15,3 in	72 cm 28,3 in	41 cm 16,1 in	23 cm 9 in	90-120 cm 35,4-47,2 in	35 kg 77,1 lb	8,8 kg 19,4 lb	8,5 kg 18,7 lb	17,3 kg 38,1 lb
<b>Medium</b>	25-36 cm 9,8-14,1 in	60-75 cm 23,6-29,5 in	46 cm 18,1 in	34-50 cm 13,3-19,6 in	30-40 cm 11,8-15,7 in	20-38 cm 7,8-14,9 in	0°-25°	40°	95-126 cm 37,4-49,6 in	55 cm 21,6 in	23-39 cm 9-15,3 in	76 cm 29,9 in	46 cm 18,1 in	28 cm 11 in	115-140 cm 45,2-55,1 in	45 kg 99,2 lb	10 kg 22 lb	12,7 kg 27,9 lb	22,7 kg 50 lb

### SEAT UNIT, the only one with a tilt-in-space system incorporated in the seat

- HEADREST incorporated in the backrest: your child's head is always supported during the tilt adjustments of the seat and of the backrest, also without accessories.
- Rigid BACKREST, reclinable and adjustable in height. The hinge coincides with the child's pelvis physiological fulcrum.
- Tilt-in-space SEAT, adjustable in width and depth. Reversible (size Small, base 869).
- Fireproof, recyclable FABRIC. Colours durable and highly resistant to sweat, light, dry cleaning and washing. Maximum resistance to rubbing. Oeko-Tex Certificate: for all fabrics that are in direct contact with the skin, to guarantee of quality and safety for health.
- Tip-up FOOTREST, flexion-extension adjustable, with legrest adjustable in height and tilt; the hinge coincides with the child's knee physiological fulcrum.

### 869 OUTDOOR four wheels base - CRASH TESTED ISO 7176/19 and ANSI RESNA WC/19

- Foldable aluminium FRAME. Push-handle adjustable in tilt: adapts to the mother's height and reduces the encumbrance of the pushchair.
- Puncture proof WHEELS, light and easy to push on every surface: front 17 cm with directional locks; rear 25 cm (size Small) or 30 cm (size Medium), with shock absorbers and independent drum brakes.

### 856 - 851 indoor high-low bases

- Foldable FRAME in steel. Adjustable in height. It is easy and safe to use, thanks to a gas spring mechanism spring, hand lever (base 851) and foot pedal (base 856)
- 10 cm swivelling WHEELS with brakes (for base 856)
- 7.5 cm swivelling WHEELS with brakes (for base 851)

ADDITIONAL COMPONENTS	<input type="checkbox"/> <b>863</b> Multiadjustable Headrest with parietal supports 	<input type="checkbox"/> <b>852</b> Headrest with parietal supports 	<b>942</b> Ergonomic headrest for Bug Small <input type="checkbox"/> S <input type="checkbox"/> M 	<b>942</b> Ergonomic headrest for Bug Medium <input type="checkbox"/> M <input type="checkbox"/> L 	<input type="checkbox"/> <b>835</b> Shaped headrest 	<input type="checkbox"/> <b>934</b> Suspension headrest attachment. Only available for 863 headrest 
	<input type="checkbox"/> <b>838</b> Adjustable lateral supports 	<input type="checkbox"/> <b>868</b> Adjustable wrappable and flexible trunk supports 	<input type="checkbox"/> <b>834</b> Padded abduction block 	<input type="checkbox"/> <b>834R</b> Adjustable padded abduction block 	<input type="checkbox"/> <b>853</b> Vest harness 	<input type="checkbox"/> <b>903</b> Five point vest harness 
	<input type="checkbox"/> <b>906</b> Five point harness 	<input type="checkbox"/> <b>894</b> 45° pelvic belt 	<input type="checkbox"/> <b>920</b> 4 Point belt 	<input type="checkbox"/> <b>824</b> Adjustable transparent tray 	<input type="checkbox"/> <b>824B</b> Wide tray just for 850- 851 -856 indoor bases. 	
	<input type="checkbox"/> <b>839</b> Adjustable front handle 	<input type="checkbox"/> <b>827</b> Foot straps 	<input type="checkbox"/> <b>892</b> Padded covering for footrest 	<input type="checkbox"/> <b>869</b> OUTDOOR four wheels base  <p>                     CRASH TESTED                      ACCORDING TO THE ISO 7176/19 and ANSI RESNA WC/19                      Bug with base 869 passed CRASH TEST according to the ISO 7176/19 and ANSI RESNA WC/19                 </p> <p>For the accessories of the 869 base, please refer to the appropriate order form</p>		