

tray and wheels included




MARKED 
 FDA registered - USA

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

F E A T U R E S

- **No edge tilting frame** in order to enable trunk rotation and side-arm motion. Smooth continuous adjustment, without snaps (through spring stopper adjuster). It's adjustable by knobs, without using wrenches, "tools free".
- **Birch wood table.** It is removable and adjustable in height, tilt and depth.
- **Trunk Support** adjustable in its circumference thanks to Velcro straps, in height, back and forward.
- **Pelvic Support** adjustable in its circumference thanks to Velcro straps, in height, back and forward.
- **Numbered bars** to find the right setting, determined by the therapist.
- **Knee-pads** adjustable in height, width, depth and in a rotatory movement.
- **Heelrests** adjustable in depth and width.
- **Birch wood footboard.** It gives a safe guide to the child. It is adjustable in height.
- **Twinned solid Wheels,** 75 mm. diam., all wheels with brakes. Thanks to the wheels the aid can be moved from one room to another having the child on it.


 **Attention!**
 Read carefully the user and maintenance manual

Please check off the interested products:



M E A S U R E S

- DONDOLINO 1 DONDOLINO 2 DONDOLINO 3




size	Max. Load	Tot.Weight	A	B	C	D	E	F	G	L	M**	N	P
1	35 kg 77,1 lb	21 kg 46,2 lb	63 cm 24,8 in	75-85 cm 29,5-33,4 in	45-79 cm 17,7-31,1 in	63 cm 24,8 in	12-26 cm 4,7-10,2 in	39-59 cm 15,3-23,2 in	51-76 cm 20-29,9 in	60 cm 23,6 in	18 cm 7 in	75-100 cm 29,5-39,3 in	0°-19°
2	45 kg 99,2 lb	24 kg 52,9 lb	77 cm 30,3 in	81-91 cm 31,8-35,8 in	52-90 cm 20,4-35,4 in	67 cm 26,3 in	12-36 cm 4,7-14,1 in	47-74 cm 18,5-29,1 in	56-90 cm 22-35,4 in	65 cm 25,5 in	23 cm 9 in	90-125 cm 35,4-49,2 in	0°-14°
3	55 kg 121,2 lb	26 kg 57,3 lb	83 cm 32,6 in	85-95 cm 33,4-37,4 in	62-107 cm 24,4-42,1 in	71 cm 27,9 in	14-48 cm 5,5-18,8 in	58-93 cm 22,8-36,6 in	70-110 cm 27,5-43,3 in	70 cm 27,5 in	28 cm 11 in	120-150 cm 47,2-59 in	0°-11°


ADDITIONAL COMPONENTS

944 Circumference reduction parts
 Thickness: 3 cm. each side



865 Multiadjustable headrest



 download and send the form to info@ormesa.com

Notes:-----

