

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



Tested by TÜV Berlin  
 FDA registered - USA  
 CE marked

**FEATURES**

- Standy** consists of a **STRONG STRUCTURE** made of tubular steel, with sliding parts chromed and rest painted with polyester and epoxy-polyester varnishes
- Removable **table** adjustable in height with semi-rigid side protection
  - **Knee pads** adjustable in height, width and depth
  - Padded **back supports** and **chest supports**
  - **Heel supports** adjustable in width and backwards-forwards
  - Non-slip aluminium **footboard**
  - Transport **wheels** allow the device to be moved from one room to another without the user in place

**MEASURES**


Please check off the interested products:  \_\_\_\_\_

**STANDY 3 medium**

Frame colour: **bronze**  
 Upholstery colour: **brown**

**STANDY 4 large**

Frame colour: **bronze**  
 Upholstery colour: **brown**



SIZE	Max load cap	Aid weight	A	B	C	D	E	F	G	H	L
<b>3 medium</b>	Kg. 90	Kg. 30	cm 89	cm 76	cm 63-93	cm 81-121	cm 69	cm 24-49	cm 74	cm 110-160	cm 30
<b>4 large</b>	Kg. 120	Kg. 38	cm 100	cm 85	cm 72-106	cm 84-129	cm 79	cm 30-57	cm 80	cm 135-180	cm 39

Attention!  
 Read carefully the user and maintenance manual

**ADDITIONAL COMPONENTS**

<input type="checkbox"/> <b>865</b> Headrest	<input type="checkbox"/> <b>857</b> Side support for the pelvis (1 piece)	<input type="checkbox"/> <b>945</b> Handle to open the back support for STANDY4	<input type="checkbox"/> <b>861</b> 4 wheel set with brake	<input type="checkbox"/> <b>902</b> Front footrests with straps
--	---	---	--	---

INGLESE - Cod 22577 - 08.2021

download and send the form to [info@ormesa.com](mailto:info@ormesa.com)