

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



Registered FDA (USA)
 Marked CE

**F
E
A
T
U
R
E
S**

- Standy Electro Wing consists of a **strong structure** made of tubular steel, painted with polyester and epoxy-polyester varnishes. Sliding parts are chromed.
- **Wheels** with brakes allow the device to be moved from one room to another with the user already in place.
- **Heel stops** adjustable in depth and abduction.
- **Knee pads** adjustable in depth, height and abduction.
- **Tray** adjustable in height, backwards-forwards and in tilt (0° - 23°).
- The **emergency button** (located under the tray) is depth adjustable using 4 screws.
- Standard **lifting band** ergonomic, with internal covering in no-slide technical fabric. Hooks to allow traction while climbing. An alternative maxi size is available (907 model).
Alternative 908 lifting band with pelvis belt. A maxi size is available (909 model).



Please check off the interested products:



**M
E
A
S
U
R
E
S**

Frame colour: **anthracite grey** Upholstery colour: **black**

- STANDY ELECTRO WING (for adults)**
- STANDARD LIFTING BAND**
- MAXI LIFTING BAND mod. 907**
- LIFTING BAND + PELVIC BELT mod. 908**
- MAXI LIFTING BAND + PELVIC BELT mod. 909**



Total Weight: **60 kg** Max. Load: **120 kg** User Height: **140-180 cm**

Attention!
 Read carefully the user and maintenance manual

ADDITIONAL COMPONENTS

866R Pelvic and backrest support and adjustable headrest



866 Pelvic, head and support and back support



867 Pelvic support



902 Front footrests with straps



INGLESE - Cod 22659 - 26.01.23