

# ORMESA®

• MADE IN ITALY SINCE 1980 •



## ORDER FORM

**BIK**  
ADAPTIVE TRICYCLE

size Extra Large

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

### F E A T U R E S

- **Ergonomic** and **modular tricycle** that helps the user to keep a proper posture and makes the ride easier thanks to its design and components.
- **Lowered frame** for easy ascent and descent.
- **Front positioned bottom bracket** for easier pedaling.
- **Pinion** to facilitate the pedaling in safety.
- **Folding frame** to save space during transport.
- **Adjustable handlebar** in height, backward-forward and reversible.
- **Wide and ergonomic handles**, to increase comfort.
- **Bell**.
- **Front and back led lamp**.
- **Front brake** with **locking device** for ascent and descent, that can be operated by the lever either on the right or on the left side of the handlebar. The distance between the brake lever and the handlebar is adjustable by screw.
- **Adjustable ergonomic saddle** in height, tilt and backward-forward
- **Rear storage basket**.
- **20" pneus wheels**. A mechanism on the rear wheels allows easy adjustment of the tension of the chain.
- **Quick release front wheel**.



Product marked **CE**




Attention!  
Read carefully the user and  
maintenance manual

Please check off the interested products:

### M E A S U R E S


#### FRAME COLOUR:

- Blue   
 Red   
 Green



#### Storage basket Internal measures

L	55,3 cm in 21,7
M	27,4 cm in 10,7
N	9,6 cm in 3,7
Max. load	35 kg lb 77



	Max. Load	Aid Weight	A	B	C1*	D	E	F	G	I	User Height (H)
XL	kg. 100 lb. 220	kg. 28 lb. 61,7	cm 64 in 25,2	cm 82 in 32,2	cm 73-87 in 28,7-34	cm 16 in 6,3	cm 168 in 66	cm 61,5 in 24,2	cm 96 in 37,7	cm 115 in 45,2	cm 155-190 in 61-74,8

\*With standard saddle. You can further adjust the height of the saddle (3cm ±) acting on the rail clamp.

Choose the configuration →

Fixed Pinion


Free Wheel

Differential and Fixed Pinion

Differential and Free Wheel

Three-speed transmission with differential. Allows braking by counter-pedaling, shifting when the tricycle is stationary, and reversing.


▼ OPTIONAL COMPONENTS, ONLY CHOOSABLE WHEN ORDERING THE TRICYCLE



- 864 Disc brake. only for models with differential and without gearbox.



- 847 Cranks adjustable in length decreases the distance between saddle and pedal of 3,5 and 5 cm.



Standard handlebar

size L     size XL

Please choose the position of the brake

Left     Right

In the absence of indications, the brake will be mounted on the right.

▼ OPTIONAL COMPONENTS, ALSO CHOOSABLE SEPARATELY

- 880 Ergonomic handles adjustable in inclination independently (1 pair). Same size as the Handlebar.



- 845 Forearm support

Left     Right  
 S-M     L-XL



- 875 Wrist retainer with support and multi-adjustable glove independently.

Left     Right  
 S-M     L-XL



- 842P Pelvic support with belt, adjustable in height.

Also select the component 883

M  
 L-XL



- 842T Trunk support, adjustable in height, inclination and backward-forward with belt and joint. Also select the component 883

M  
 L-XL



- 883 Rod adjustable in height for Thoracic/Pelvic support (single rod). Please select even the type of support.



- 877 Flat backrest



- 843 Abductive saddle with holder for supports



- 876 Monocycle saddle with holder for supports



- 840 Pedal cages with straps (1 pair).



- 841 Pedal cage with length and width adjustable foot imprints and fastening straps, with support kit.




- 846 Ankle support adjustable in height, with support kit. Only with component 841.



- 881 Steering limiter.



- 882 Steering limiter with springs to return the steering to the straight position.









- 879 Front basket.



- 878 878 MTB tires (3 pcs.) in alternative to standard tyres



COMPONENTS MEASURES	842P	842T																											
																													
	<table border="1"> <tr><th></th><th>A</th></tr> <tr><td>mini/s</td><td>8,5-14 cm/in 3,3-5,5</td></tr> <tr><td>M</td><td>8,5-21 cm/in 3,3-8,3</td></tr> <tr><td>L/XL</td><td>8,5-23 cm/in 3,3-9</td></tr> </table>		A	mini/s	8,5-14 cm/in 3,3-5,5	M	8,5-21 cm/in 3,3-8,3	L/XL	8,5-23 cm/in 3,3-9	<table border="1"> <tr><th></th><th>A</th></tr> <tr><td>mini/s</td><td>19-30,5 cm/in 7,4-12</td></tr> <tr><td>M</td><td>19-40 cm/in 7,4-15,7</td></tr> <tr><td>L/XL</td><td>19-46 cm/in 7,4-18,1</td></tr> </table>		A	mini/s	19-30,5 cm/in 7,4-12	M	19-40 cm/in 7,4-15,7	L/XL	19-46 cm/in 7,4-18,1											
	A																												
mini/s	8,5-14 cm/in 3,3-5,5																												
M	8,5-21 cm/in 3,3-8,3																												
L/XL	8,5-23 cm/in 3,3-9																												
	A																												
mini/s	19-30,5 cm/in 7,4-12																												
M	19-40 cm/in 7,4-15,7																												
L/XL	19-46 cm/in 7,4-18,1																												
<b>842P - 842T</b>	<b>843</b>	<b>876</b>																											
																													
<table border="1"> <tr><th></th><th>B</th><th>C</th></tr> <tr><td>mini/s</td><td>10,5 cm/in 4,1</td><td>23 cm/in 9</td></tr> <tr><td>M</td><td>12,5 cm/in 4,9</td><td>31 cm/in 12,2</td></tr> <tr><td>L/XL</td><td>17,5 cm/in 6,8</td><td>35 cm/in 13,7</td></tr> </table>		B	C	mini/s	10,5 cm/in 4,1	23 cm/in 9	M	12,5 cm/in 4,9	31 cm/in 12,2	L/XL	17,5 cm/in 6,8	35 cm/in 13,7	<table border="1"> <tr><th></th><th>A</th><th>B</th></tr> <tr><td>mini/s</td><td>31 cm/in 12,2</td><td>12,5 cm/in 4,9</td></tr> <tr><td>M/L/XL</td><td>40 cm/in 15,7</td><td>14 cm/in 5,5</td></tr> </table>		A	B	mini/s	31 cm/in 12,2	12,5 cm/in 4,9	M/L/XL	40 cm/in 15,7	14 cm/in 5,5	<table border="1"> <tr><th>A</th><td>31 cm 12,2"</td></tr> <tr><th>B</th><td>13,5 cm 5,3"</td></tr> <tr><th>C</th><td>9 cm 3,5"</td></tr> </table>	A	31 cm 12,2"	B	13,5 cm 5,3"	C	9 cm 3,5"
	B	C																											
mini/s	10,5 cm/in 4,1	23 cm/in 9																											
M	12,5 cm/in 4,9	31 cm/in 12,2																											
L/XL	17,5 cm/in 6,8	35 cm/in 13,7																											
	A	B																											
mini/s	31 cm/in 12,2	12,5 cm/in 4,9																											
M/L/XL	40 cm/in 15,7	14 cm/in 5,5																											
A	31 cm 12,2"																												
B	13,5 cm 5,3"																												
C	9 cm 3,5"																												
<b>841</b>	<b>Standard handlebar</b>																												
																													
<table border="1"> <tr><th colspan="3">841 foot imprint size</th></tr> <tr><th>A</th><th>B</th><th>C</th></tr> <tr><td>mini/s</td><td>from 5 to 7 cm from 1,9 to 2,7 in</td><td>from 12,5 to 16,5 cm from 4,9 to 6,5 in</td><td>from 6 to 7 cm from 2,3 to 2,7 in</td></tr> <tr><td>M/L/XL</td><td>from 8,2 to 9,4 cm from 3,2 to 3,7 in</td><td>from 20 to 22,5 cm from 7,8 to 8,8 in</td><td>from 10,4 to 12 cm from 4 to 4,7 in</td></tr> </table>	841 foot imprint size			A	B	C	mini/s	from 5 to 7 cm from 1,9 to 2,7 in	from 12,5 to 16,5 cm from 4,9 to 6,5 in	from 6 to 7 cm from 2,3 to 2,7 in	M/L/XL	from 8,2 to 9,4 cm from 3,2 to 3,7 in	from 20 to 22,5 cm from 7,8 to 8,8 in	from 10,4 to 12 cm from 4 to 4,7 in	<table border="1"> <tr><th>A</th></tr> <tr><td>mini</td><td>45 cm / in 17,7</td></tr> <tr><td>S</td><td>50 cm / in 19,6</td></tr> <tr><td>M</td><td>55 cm / in 21,6</td></tr> <tr><td>L</td><td>60 cm / in 23,6</td></tr> <tr><td>XL</td><td>64 cm / in 25,2</td></tr> </table>	A	mini	45 cm / in 17,7	S	50 cm / in 19,6	M	55 cm / in 21,6	L	60 cm / in 23,6	XL	64 cm / in 25,2			
841 foot imprint size																													
A	B	C																											
mini/s	from 5 to 7 cm from 1,9 to 2,7 in	from 12,5 to 16,5 cm from 4,9 to 6,5 in	from 6 to 7 cm from 2,3 to 2,7 in																										
M/L/XL	from 8,2 to 9,4 cm from 3,2 to 3,7 in	from 20 to 22,5 cm from 7,8 to 8,8 in	from 10,4 to 12 cm from 4 to 4,7 in																										
A																													
mini	45 cm / in 17,7																												
S	50 cm / in 19,6																												
M	55 cm / in 21,6																												
L	60 cm / in 23,6																												
XL	64 cm / in 25,2																												

Notes:-----

-----

-----

-----

-----

-----

-----

-----

-----