



MARKED

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

FEATURES

- **Frame** height-adjustable, made of steel painted with non-toxic epoxy paint.
- Fire-retardant **upholstery**.
- **Seat** triangular to hold the pelvis (seat/back at a 90° angle). It is padded and wrap-around, upholstered with strong, washable plastic.
- **Bumpers** elastic wall-protector bumpers and protective balls on the upper ends of the rods.
- **Wheels** 80 mm fully-swivelling, rolling on ball-bearings, 2 with brakes.
- **Handle bar** adjustable in both height and depth.
- **Abduction block** adjustable to contain the pelvis.



Attention!
Read carefully the user and maintenance manual

Please check off the interested products



BIRILLO 1 Extra Small

Frame colour: **blue**


Upholstery colours: **yellow-green-red**

BIRILLO 2 Small

BIRILLO 3 Medium

BIRILLO 4 Large

MEASURES



SIZE	Aid Weight	Turning Width	Max. Cap. Load	A	B	C	D	E	F	G	H
1	8,5 kg 18,7 lb	67 cm 26,3 in	30 kg 66,1 lb	18 cm 7 in	59 cm 23,2 in	50 cm 19,6 in	20-28 cm 7,8-11 in	34-48 cm 13,3-18,8 in	18 cm 7 in	22 cm 8,6 in	70-95 cm 27,5-37,4 in
2	10 kg 22 lb	73 cm 28,7 in	40 kg 88,1 lb	19 cm 7,4 in	65 cm 25,5 in	51 cm 20 in	23-31 cm 9-12,2 in	40-56 cm 15,7-22 in	20 cm 7,8 in	26 cm 10,2 in	85-105 cm 33,4-41,3 in
3	11 kg 24,2 lb	81 cm 31,8 in	50 kg 110,2 lb	21 cm 8,2 in	73 cm 28,7 in	59 cm 23,2 in	23-35 cm 9-13,7 in	40-65 cm 15,7-25,5 in	23 cm 9 in	30 cm 11,8 in	95-120 cm 37,4-47,2 in
4	12,8 kg 28,2 lb	89 cm 35 in	60 kg 132,2 lb	26 cm 10,2 in	78 cm 30,7 in	60 cm 23,6 in	22-40 cm 8,6-15,7 in	43-74 cm 16,9-29,1 in	30 cm 11,8 in	38 cm 14,9 in	110-140 cm 43,3-55,1 in

ADDITIONAL COMPONENTS

<input type="checkbox"/> 812 Directional lock 	<input type="checkbox"/> 824 Transparent tray 	<input type="checkbox"/> 86 Headrest
---	---	--

Notes:

INGLESE - Cod 22576 -28.10.24