

• MADE IN ITALY SINCE 1980 •



## ORDER FORM

**BIK**  
ADAPTIVE TRICYCLE

size Extra Large

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

### F E A T U R E S

- **Ergonomic** and **modular tricycle** that helps the user to keep a proper posture and makes the ride easier thanks to its design and components.
- **Lowered frame** for easy ascent and descent.
- **Front positioned bottom bracket** for easier pedaling.
- **Pinion** to facilitate the pedaling in safety.
- **Folding frame** to save space during transport.
- **Adjustable handlebar** in height, backward-forward and reversible.
- **Wide and ergonomic handles**, to increase comfort.
- **Bell**.
- **Front and back led lamp**.
- **Front brake** with **locking device** for ascent and descent, that can be operated by the lever either on the right or on the left side of the handlebar. The distance between the brake lever and the handlebar is adjustable by screw.
- **Adjustable ergonomic saddle** in height, tilt and backward-forward
- **Rear storage basket**.
- **20" pneus wheels**. A mechanism on the rear wheels allows easy adjustment of the tension of the chain.
- **Quick release front wheel**.



Product marked




Attention!  
Read carefully the user and  
maintenance manual

Please check off the interested products:

### M E A S U R E S


#### FRAME COLOUR:

- Blue
- Red
- Green



#### Storage basket Internal measures

L	55,3 cm in 21,7
M	27,4 cm in 10,7
N	9,6 cm in 3,7
Max. load	35 kg lb 77



	Max. Load	Aid Weight	A	B	C1*	D	E	F	G	I	User Height (H)
XL	kg. 100 lb. 220	kg. 28 lb. 61,7	cm 64 in 25,2	cm 82 in 32,2	cm 73-87 in 28,7-34	cm 16 in 6,3	cm 168 in 66	cm 61,5 in 24,2	cm 96 in 37,7	cm 115 in 45,2	cm 155-190 in 61-74,8

\*With standard saddle. You can further adjust the height of the saddle (3cm ±) acting on the rail clamp.





### Choose the configuration


<input type="checkbox"/> Fixed Pinion	<input type="checkbox"/> Free Wheel	<input type="checkbox"/> Differential and Fixed Pinion	<input type="checkbox"/> Differential and Free Wheel	<input type="checkbox"/> Differential and "Switch" system. It allows you to switch from freewheel to fixed pinion and vice versa.	<input type="checkbox"/> Three-speed transmission with differential. Allows braking by counter-pedaling, shifting when the tricycle is stationary, and reversing.
---------------------------------------	-------------------------------------	--	--	---	---

OPTIONAL  
COMPONENTS,  
ONLY  
CHOOSABLE  
WHEN ORDERING  
THE TRICYCLE

<input type="checkbox"/> 864 Disc brake. only for models with differential and without gearbox. Not available with "Switch" system.	Standard handlebar <input type="checkbox"/> size L <input type="checkbox"/> size XL <input type="checkbox"/> Left <input type="checkbox"/> Right	<input type="checkbox"/> 847 Cranks adjustable in length decreases the distance between saddle and pedal of 3,5 and 5 cm.
Please choose the position of the brake. In the absence of indications, the brake will be mounted on the right.		

### OPTIONAL COMPONENTS, ALSO CHOOSABLE SEPARATELY

<input type="checkbox"/> 880 Ergonomic handles adjustable in inclination independently (1 pair). Same size as the Handlebar.  	<input type="checkbox"/> 845 Forearm support. Size L-XL.    <input type="checkbox"/> Left <input type="checkbox"/> Right	<input type="checkbox"/> 875 Wrist retainer with support and multi-adjustable glove independently. Size L-XL.    <input type="checkbox"/> Left <input type="checkbox"/> Right	<input type="checkbox"/> 842P Pelvic support with belt, adjustable in height. Also select the component 883.    <input type="checkbox"/> M <input type="checkbox"/> L-XL
<input type="checkbox"/> 842T Trunk support, adjustable in height, inclination and backward-forward with belt and joint. Also select the component 883.    <input type="checkbox"/> M <input type="checkbox"/> L-XL	<input type="checkbox"/> 883 Rod adjustable in height for Thoracic/Pelvic support (single rod). Please select even the type of support.  	<input type="checkbox"/> 877 Flat backrest.  	<input type="checkbox"/> 843 Abductive saddle with holder for supports.  
<input type="checkbox"/> 876 Monocycle saddle with holder for supports.  	<input type="checkbox"/> 103 Monocycle saddle padding.  	<input type="checkbox"/> 112 Comfortable larger saddle, with wider support surface.  	
<input type="checkbox"/> 840 Pedal cages with straps (1 pair).  	<input type="checkbox"/> 841 Pedal cage with lenght and width adjustable foot imprints and fastening straps, with support kit.    <input type="checkbox"/> Left <input type="checkbox"/> Right	<input type="checkbox"/> 846 Ankle support adjustable in height, with support kit. Only with component 841.    <input type="checkbox"/> Left <input type="checkbox"/> Right  <input type="checkbox"/> S-M <input type="checkbox"/> L-XL	



COMPONENTS MEASURES	842P	842T	842P - 842T																											
	<table border="1"> <tr> <td>mini/s</td> <td>8,5-14 cm / 3,3-5,5 in</td> </tr> <tr> <td>M</td> <td>8,5-21 cm / 3,3-8,3 in</td> </tr> <tr> <td>L/XL</td> <td>8,5-23 cm / 3,3-9 in</td> </tr> </table>	mini/s	8,5-14 cm / 3,3-5,5 in	M	8,5-21 cm / 3,3-8,3 in	L/XL	8,5-23 cm / 3,3-9 in	<table border="1"> <tr> <td>mini/s</td> <td>19-30,5 cm / 7,4-12 in</td> </tr> <tr> <td>M</td> <td>19-40 cm / 7,4-15,7 in</td> </tr> <tr> <td>L/XL</td> <td>19-46 cm / 7,4-18,1 in</td> </tr> </table>	mini/s	19-30,5 cm / 7,4-12 in	M	19-40 cm / 7,4-15,7 in	L/XL	19-46 cm / 7,4-18,1 in	<table border="1"> <tr> <td>mini/s</td> <td>10,5 cm / 4,1 in</td> <td>23 cm / 9 in</td> </tr> <tr> <td>M</td> <td>12,5 cm / 4,9 in</td> <td>31 cm / 12,2 in</td> </tr> <tr> <td>L/XL</td> <td>17,5 cm / 6,8 in</td> <td>35 cm / 13,7 in</td> </tr> </table>	mini/s	10,5 cm / 4,1 in	23 cm / 9 in	M	12,5 cm / 4,9 in	31 cm / 12,2 in	L/XL	17,5 cm / 6,8 in	35 cm / 13,7 in						
mini/s	8,5-14 cm / 3,3-5,5 in																													
M	8,5-21 cm / 3,3-8,3 in																													
L/XL	8,5-23 cm / 3,3-9 in																													
mini/s	19-30,5 cm / 7,4-12 in																													
M	19-40 cm / 7,4-15,7 in																													
L/XL	19-46 cm / 7,4-18,1 in																													
mini/s	10,5 cm / 4,1 in	23 cm / 9 in																												
M	12,5 cm / 4,9 in	31 cm / 12,2 in																												
L/XL	17,5 cm / 6,8 in	35 cm / 13,7 in																												
843	<table border="1"> <tr> <td>A</td> <td>B</td> </tr> <tr> <td>mini/s</td> <td>31 cm / 12,2 in</td> <td>12,5 cm / 4,9 in</td> </tr> <tr> <td>M/L/XL</td> <td>40 cm / 15,7 in</td> <td>14 cm / 5,5 in</td> </tr> </table>	A	B	mini/s	31 cm / 12,2 in	12,5 cm / 4,9 in	M/L/XL	40 cm / 15,7 in	14 cm / 5,5 in	<table border="1"> <tr> <td>A</td> <td>C</td> <td>B</td> </tr> <tr> <td>mini/s</td> <td>31 cm / 12,2 in</td> <td>13,5 cm / 5,3 in</td> <td>9 cm / 3,5 in</td> </tr> <tr> <td>M/L/XL</td> <td>40 cm / 15,7 in</td> <td>14 cm / 5,5 in</td> <td>9 cm / 3,5 in</td> </tr> </table>	A	C	B	mini/s	31 cm / 12,2 in	13,5 cm / 5,3 in	9 cm / 3,5 in	M/L/XL	40 cm / 15,7 in	14 cm / 5,5 in	9 cm / 3,5 in	<table border="1"> <tr> <td>A</td> <td>B</td> </tr> <tr> <td>mini/s</td> <td>25 cm / 9,8 in</td> <td>17 cm / 6,7 in</td> </tr> <tr> <td>M/L/XL</td> <td>25,5 cm / 10 in</td> <td>20 cm / 7,8 in</td> </tr> </table>	A	B	mini/s	25 cm / 9,8 in	17 cm / 6,7 in	M/L/XL	25,5 cm / 10 in	20 cm / 7,8 in
A	B																													
mini/s	31 cm / 12,2 in	12,5 cm / 4,9 in																												
M/L/XL	40 cm / 15,7 in	14 cm / 5,5 in																												
A	C	B																												
mini/s	31 cm / 12,2 in	13,5 cm / 5,3 in	9 cm / 3,5 in																											
M/L/XL	40 cm / 15,7 in	14 cm / 5,5 in	9 cm / 3,5 in																											
A	B																													
mini/s	25 cm / 9,8 in	17 cm / 6,7 in																												
M/L/XL	25,5 cm / 10 in	20 cm / 7,8 in																												
841	<table border="1"> <tr> <td colspan="3">841 foot imprint size</td> </tr> <tr> <td>A</td> <td>B</td> <td>C</td> </tr> <tr> <td>mini/s</td> <td>from 5 to 7 cm from 1,9 to 2,7 in</td> <td>from 12,5 to 16,5 cm from 4,9 to 6,5 in</td> <td>from 6 to 7 cm from 2,3 to 2,7 in</td> </tr> <tr> <td>M/L/XL</td> <td>from 8,2 to 9,4 cm from 3,2 to 3,7 in</td> <td>from 20 to 22,5 cm from 7,8 to 8,8 in</td> <td>from 10,4 to 12 cm from 4 to 4,7 in</td> </tr> </table>	841 foot imprint size			A	B	C	mini/s	from 5 to 7 cm from 1,9 to 2,7 in	from 12,5 to 16,5 cm from 4,9 to 6,5 in	from 6 to 7 cm from 2,3 to 2,7 in	M/L/XL	from 8,2 to 9,4 cm from 3,2 to 3,7 in	from 20 to 22,5 cm from 7,8 to 8,8 in	from 10,4 to 12 cm from 4 to 4,7 in	<table border="1"> <tr> <td>A</td> </tr> <tr> <td>mini</td> <td>45 cm / 17,7 in</td> </tr> <tr> <td>S</td> <td>50 cm / 19,6 in</td> </tr> <tr> <td>M</td> <td>55 cm / 21,6 in</td> </tr> <tr> <td>L</td> <td>60 cm / 23,6 in</td> </tr> <tr> <td>XL</td> <td>64 cm / 25,2 in</td> </tr> </table>	A	mini	45 cm / 17,7 in	S	50 cm / 19,6 in	M	55 cm / 21,6 in	L	60 cm / 23,6 in	XL	64 cm / 25,2 in			
841 foot imprint size																														
A	B	C																												
mini/s	from 5 to 7 cm from 1,9 to 2,7 in	from 12,5 to 16,5 cm from 4,9 to 6,5 in	from 6 to 7 cm from 2,3 to 2,7 in																											
M/L/XL	from 8,2 to 9,4 cm from 3,2 to 3,7 in	from 20 to 22,5 cm from 7,8 to 8,8 in	from 10,4 to 12 cm from 4 to 4,7 in																											
A																														
mini	45 cm / 17,7 in																													
S	50 cm / 19,6 in																													
M	55 cm / 21,6 in																													
L	60 cm / 23,6 in																													
XL	64 cm / 25,2 in																													

Notes: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_