

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



Product marked **CE**

FEATURES

- **Ergonomic and modular tricycle** that helps the user to keep a proper posture and makes the ride easier thanks to its design and components.
- **Lowered frame** for easy ascent and descent.
- **Front positioned bottom bracket** for easier pedaling.
- **Optimized pinion** to facilitate the pedaling in safety.
- **Adjustable handlebar** in height, backward-forward and reversible.
- **Wide and ergonomic handles**, to increase comfort.
- **Bell**.
- **Back led lamp**.
- **Cranks adjustable** in lenght.
- **Front brake with locking device** for ascent and descent, that can be operated by the thumb-lock stopper either on the right or on the left side of the handlebar. The distance between the brake lever and the handlebar is adjustable by screw.
- **Adjustable ergonomic saddle** in height, tilt and backward-forward
- **Rear storage basket**.
- **12" solid wheels**. A mechanism on the rear wheels allows easy adjustment of the tension of the chain.



Attention!
Read carefully the user and maintenance manual


Please check off the interested products: ☒

MEASURES

FRAME COLOUR:


- ☐ Blue
- ☐ Red
- ☐ Pink

View from above




Storage basket
Internal measures

L	41 cm	In 16
M	18 cm	In 7
N	6,6 cm	In 2,6
Max. Load	15 kg	Lb 33



$$C1 = C2 - 2cm$$

NEW



C2= leg length

	Max. load	Aid weight	A	B	C1*	D	E	F	Use height(H)
size mini	Kg. 30 lb. 66,1	Kg. 15,5 lb. 34,1	cm 45 in. 17,7	cm 62 in. 24,4	cm 39-53 in. 15,3-20,8	cm 12 in. 4,7	cm 103 in. 40,5	cm 36 in. 14,2	cm 90-110 in. 35,4-43,3

*With standard saddle. You can further adjust the height of the saddle (3cm ±) acting on the clamp.

CHOOSE THE CONFIGURATION >



Fixed Pinion



Free Wheel



"Switch" system.

It allows you to switch from freewheel to fixed pinion and vice versa.



Please **CHOOSE THE POSITION OF THE BRAKE** on Standard Handlebar or on Ergonomic Handles 880.

☐ Left

☐ Right

In the absence of indications, the brake will be mounted on the right.

OPTIONAL COMPONENTS, ALSO CHOOSABLE SEPARATELY

☐ **880 Ergonomic handles adjustable in inclination independently (1 pair).** Same size as the Handlebar.

Choose the position of the brake on Ergonomic Handles



845 Forearm support.

☐ Left

☐ Right



875 Wrist retainer with support and multi-adjustable glove independently.

☐ Left

☐ Right



☐ **842P Pelvic support with belt, adjustable in height.**

Also select the component 883.



☐ **842T Trunk support, adjustable in height, inclination and backward-forward with belt and joint.** Also select the component 883.



☐ **883 Rod adjustable in height for Thoracic/Pelvic support (single rod).** Please select even the type of support.



☐ **877 Flat backrest.**



☐ **843 Abductive saddle with holder for supports.**



☐ **876 Monocycle saddle with holder for supports.**



☐ **103 Monocycle saddle padding.**



☐ **112 Larger saddle, with wider support surface.**



NEW


☐ **840 Pedal cages with straps (1 pair).**



841 Pedal cage with length and width adjustable foot imprints and fastening straps, with support kit.

☐ Left


☐ Right



846 Ankle support adjustable in height, with support kit. Only with component 841.

☐ Left

☐ Right



☐ **923 Additional assistant handle.**



☐ **963** Steering assistant handle adjustable in height. Cannot be installed with 881 or 882 components.



☐ **881** Steering limiter.



☐ **882** Steering limiter with springs to return the steering to the straight position.




☐ **879** Front basket.




COMPONENTS MEASURES

842P




A	
mini/s	8,5-14 cm / 3,3-5,5 in
M	8,5-21 cm / 3,3-8,3 in
L/XL	8,5-23 cm / 3,3-9 in

842T




A	
mini/s	19-30,5 cm / 7,4-12 in
M	19-40 cm / 7,4-15,7 in
L/XL	19-46 cm / 7,4-18,1 in

842P - 842T




	B	C
mini/s	10,5 cm / 4,1 in	23 cm / 9 in
M	12,5 cm / 4,9 in	31 cm / 12,2 in
L/XL	17,5 cm / 6,8 in	35 cm / 13,7 in

843




	A	B
mini/s	31 cm / 12,2 in	12,5 cm / 4,9 in
M/L/XL	40 cm / 15,7 in	14 cm / 5,5 in

876




A	31 cm / 12,2 in
B	13,5 cm / 5,3 in
C	9 cm / 3,5 in

Standard saddle




	A	B
mini/s	25 cm / 9,8 in	17 cm / 6,7 in
M/L/XL	25,5 cm / 10 in	20 cm / 7,8 in

112




	A	B
mini/s	25,5 cm / 10 in	20 cm / 7,8 in
M/L/XL	25,5 cm / 10 in	24,5 cm / 9,6 in

841



841 foot imprint size			
	A	B	C
mini/s	from 5 to 7 cm from 1,9 to 2,7 in	from 12,5 to 16,5 cm from 4,9 to 6,5 in	from 6 to 7 cm from 2,3 to 2,7 in
M/L/XL	from 8,2 to 9,4 cm from 3,2 to 3,7 in	from 20 to 22,5 cm from 7,8 to 8,8 in	from 10,4 to 12 cm from 4 to 4,7 in

Standard handlebar



	A
mini	45 cm / 17,7 in
S	50 cm / 19,6 in
M	55 cm / 21,6 in
L	60 cm / 23,6 in
XL	64 cm / 25,2 in

Notes: _____