

DISTRIBUTOR/AGENT	DATE OF ORDER
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



ISO 7176/19 and  
ANSI RESNA WC/19



very breathable  
and fire retardant  
upholstery



Marked



Attention!  
Read carefully the user and  
maintenance manual

## FEATURES

- **More sizes in one frame**
- **Frame:** light and umbrella folding system, in aluminium. It is painted with non-toxic epoxy powders. When folded Trollino is compact and it allows you to carry it as if it was a travel trolley.
- **Push handles** with **manual and parking brakes**.
- Resistant and durable **backrest** guaranteed by the cross on the back. Reclinable on three positions and easily adjustable in height.
- Slightly tilted **seat**. It is easily adjustable in depth and width thanks to padded pelvic side supports (additional components 958). The cross below ensures stability, comfort and durability.
- Extremely breathable **upholstery** (not removable). It is resistant and easily washable. **Fire Retardant** according to *EN 1021-1:2014* and *EN 1021-2:2014*.
- User manual **bag** and wrench on the back of the backrest.
- **Legrest** adjustable in height.
- Elevatable and swivel **footrests**, without encumbrance.
- **Solid wheels:** 17 cm (7 in) front swiveling wheels and 25 cm (10 in) rear wheels with **independent drum brakes**. **Overcoming step bar** on both rear wheels.

## Please check off the interested products



## MEASURES


Backrest colour:      Seat colour: **Black**      Frame colour: **Black**

- ☐ Red  
☐ Pink  
☐ Blue  
☐ Black

Unit of Measurement: cm (in)



Total Weight: 14,8 kg (32,6 lb)  
Max. Load: 55 kg (121 lb)



**A:** 36 cm (14 in). For different width (32 and 28 cm - 12,6 and 11 in) please select 958 padded pelvic side supports.  
\*Height adjustment range with add-on 951 max. 105 cm (41,3 in).

☐ **951** Height adjustable push handles. Tools free!  
COMPONENT TO BE ORDERED WITH THE PRODUCT



### COMPONENTS TO CHOOSE ALSO AFTERWARDS

**835** Shaped and  
enveloping headrest  
adjustable in height.

☐ S-M




**865** Padded  
headrest adjustable in  
height.

☐ S-M




**838** Trunk side  
supports adjustable in  
height, width and  
inclination.

☐ mini ☐ S



**868** Wrappable and  
flexible trunk support,  
multiadjustable in height,  
width and inclination.

☐ S



☐ **958** Padded  
pelvic side supports  
32 cm (12,6 in) seat width.

☐ **Extra Padding**  
28 cm (11 in) seat width




**947** Pelvic belt with  
variable angle.

☐ S ☐ M




**920** 4 point pelvic.  
belt

☐ S ☐ M




**853** Vest harness.

☐ mini ☐ S




**853 slim**  
4 point shaped  
harness.


☐ M



☐ **834N** Padded  
and removable  
abduction block.



☐ **839**  
Front handlebar  
adjustable in tilt.



☐ **960**  
Heel rests.



☐ **956** Footrests  
adjustable in tilt.



☐ **858**  
Shopping basket.



☐ **819** Elevating  
Canopy with bag.  
Foldable with the  
frame.



☐ **825** Rain cover  
to fix on the 819 canopy.



☐ **825B** Back  
and lateral protection  
clothe to fix on the 819  
canopy.




☐ **818** Thermic  
cover.



☐ **968** Padded  
covers for armrests.  
Easily removable for  
washing.



☐ **844** Padded  
side protection  
covers.












☐ **812** Front  
wheel direction  
locks.



☐ **891** Set of 4 tie down hooks  
for the forward-facing transport of the pushchair in a  
motor vehicle.



COMPONENTS MEASURES	<b>947</b>  <p>small A=19 cm/7,4 in B=6 cm/2,3 in</p> <p>medium A=25,5 cm/10 in B=8 cm/3,1 in</p>	<b>920</b>  <p>small A=16 cm/6,2 in B=15,5 cm/6,1 in</p> <p>medium A=23 cm/9 in B=19 cm/7,4 in</p>	<b>838</b>  <p>mini A=9 cm/3,5 in B=7 cm/2,7 in</p> <p>small A=10 cm/3,9 in B=12 cm/4,7 in</p>	 <p>C= 14-23 cm/5,5-9 in</p>
	<b>868</b>  <p>small A=26 cm/10,2 in B=9 cm/3,5 in</p>	<b>853</b>  <p>mini A=20,5 cm/8 in B=28 cm/11 in</p> <p>small A=25 cm/9,8 in B=34 cm/13,3 in</p>	<b>853 slim</b>  <p>medium A=33 cm/12,9 in B=46 cm/18,1 in C=5 cm/1,9 in</p>	<b>834N</b>  <p>mini/S A=12 cm/4,7 in B=4 cm/1,5 in C=14 cm/5,5 in</p>
	<b>951</b>  <p>7 cm 2,7 in</p>			

Notes.....

.....

.....

.....

.....

.....

.....

.....

.....

.....

.....