

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

FEATURES

Standy consists of a **STRONG STRUCTURE** made of tubular steel, with sliding parts chromed and rest painted with polyester and epoxi-polyester varnishes

- Removable **table** adjustable in height with semi-rigid side protection
- **Knee pads** adjustable in height, width and depth
- Padded **back supports** and **chest supports**
- **Heel supports** adjustable in width and backwards-forwards
- Non-slip aluminium **footboard**
- Transport **wheels** allow the device to be moved from one room to another without the user in place



STANDY 4
large



STANDY 3
medium

Tested by TÜV Berlin
FDA registered - USA
CE marked


Please check off the interested products:



Frame colour: **bronze**

Upholstery colour: **brown**


☐ STANDY 3



Attention!
Read carefully the user
and maintenance
manual

SIZE	Max. Load	Tot. Weight	A	B	C	D	E	F	G	H	L
STANDY 3	90 kg 198,4 lb	30 kg 66,1 lb	89 cm 35 in	76 cm 29,2 in	63-93 cm 24,8-36,6 in	81-121 cm 31,8-47,6 in	69 cm 27,1 in	24-49 cm 9,4-19,2 in	74 cm 29,1 in	110-160 cm 43,3-62,9 in	30 cm 11,8 in
STANDY 4	120 kg 264,5 lb	38 kg 83,7 lb	100 cm 39,3 in	85 cm 33,4 in	72-106 cm 28,3-41,7 in	84-129 cm 33-50,7 in	79 cm 31,1 in	30-57 cm 11,8-22,4 in	80 cm 31,4 in	135-180 cm 53,1-70,8 in	39 cm 15,3 in

☐ STANDY 4



865 Headrest



857 Side support for the pelvis

 1 pc.

 2 pc.




 **945** Handle to open the back support, for Standy 4



861 4 wheel set with brake



902 Front footrests with straps



Notes: