

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

FUNCTIONS

- Light and foldable **frame** for transport and storage.
- Ergonomic **table top**, with large surface and recess. Made of a highly resistant material, easy to sanitize and hygienic.
- Simple, reliable and precise **height adjustment system**, with internal **friction system** that prevents the accidental fall of the table top.
- **Tools free**, adjustments without the use of tools.
- Provided of **bookend**, **pen holder** and **lateral hook**.
- It is available in **three sizes** for children, teenagers and adults.



Marked **CE**

Please check off the interested products:



- ☐ Multifunction Table **size Small**
- ☐ Multifunction Table **size Medium**
- ☐ Multifunction Table **size Large**


Frame colour: **Pebble grey**

Table top colour: **Silk grey**



Attention!
Read carefully the user and
maintenance manual

MEASURES



	A	B	C	D	E	F	G	I	L	User height	Total weight
Small	64 cm 25,1 in	63 cm 24,8 in	54 cm 21,2 in	from 44 to 57 cm from 17,3 to 22,4 in	17 cm 6,6 in	84 cm 33 in	from 0° to 60°	38,3 cm 15 in	16 cm 6,2 in	from 75 to 90 cm from 29,5 to 35,4 in	6,7 kg 14,7 lb
Medium	80 cm 31,4 in	80 cm 31,4 in	61 cm 24 in	from 53 to 77 cm from 20,8 to 30,3 in	16 cm 6,2 in	104 cm 40,9 in	from 0° to 60°	42,4 cm 16,6 in	14 cm 5,5 in	from 85 to 120 cm from 33,4 to 47,2 in	8,2 kg 18 lb
Large	80 cm 31,4 in	80 cm 31,4 in	61 cm 24 in	from 69 to 108 cm from 27,1 to 42,5 in	16 cm 6,2 in	120 cm 47,2 in	from 0° to 60°	42,4 cm 16,6 in	14 cm 5,5 in	from 115 to 170 cm from 45,2 to 66,9 in	9,2 kg 20,2 lb

Notes: