

Seating system with double cross frame

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

FUNCTIONS

- Lightweight aluminium **frame**, painted with non-toxic epoxy powders. It folds just like an umbrella.
- **Semi-rigid backrest**, reclines in three positions with zippers.
- **Height adjustable Handles, Tools free.**
- **Semi-rigid seat bottom**, slightly inclined.
- **Padded upholstery** removable for washing with **removable aluminium reinforcement ribs**. Extremely breathable and quick to dry, excellent for **load distribution**; **Fireproof fabric** according to EN 1021-1:2014 and EN 1021-2:2014.
- **Net basket** on the back of the backrest.
- Removable **footrest** height adjustable and foldable.
- **Legrest strap.**
- Puncture-proof **wheels**: front 17,5 cm swivelling, rear 25 cm with independent brakes **and overcoming step bar.**



CE marked





Attention!
Read carefully the user and maintenance manual

Please check off the interested products: ☒ _____

Upholstery colour: **Black and grey** Frame colour: **Black**

☐ **CLIP 4**

Measures in cm (Inches):



Max. load	Total weight	A	B*	C	D	E	F	G	H	I	L	M	N	O	P	Q	R
76 kg 167,5 lb	13,9 kg 30,6 lb	45 (17,7)	45 (17,7)	57 (22,4)	113 (44,4)	0°-10°	38 (14,9)	71 (27,9)	20°	59 (23,2)	37-47 (14,5-18,5)	89 (35)	108 (42,5)	117 (46)	60 (23,6)	37 (14,5)	150-170 (59-66,9)

* **Seat width** without padded sides (B)

☐ **835** Headrest with parietal supports.



☐ **816** Trunk lateral containment pads.



☐ **853** Vest harness. (size M)



☐ **906** Five point harness. (size M)



☐ **903** Five point vest harness. (size M)



☐ **894** 45° Pelvic belt. (size M)



☐ **834N** Padded and removable abduction block. (size M)



☐ **858** Shopping basket.



☐ **827** Foot straps.



☐ **812** Front wheel direction locks.



☐ **819** Canopy.




☐ **818** Thermal cover.



☐ **825** Rain cover. (to fix to the canopy 819).




853



Medium
A=39,5 cm/15,5 in
B=39 cm/15,3 in

834N



Medium
A=15 cm/5,9 in
B=4 cm/1,5 in
C=17 cm/6,6 in

Notes.....
.....
.....
.....