

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

FEATURES

- **Frame** painted with non-toxic epoxy powders. Compact and adjustable in height - *to allow access to different table tops and work surfaces* - thanks to the gas spring for the adjustment with the child seated. The **international height adjustment patented system**, typical of the "GRILLO product family", allows for a compact, space-saving base while offering a wide range of height adjustments.
- Grillo Adaptive Seating has a very **breathable, comfortable, durable, easy to wash and quick to dry** fabric, composed by:
 - » Load-bearing upholstery, not removable: the essential design texture allows excellent breathability that **reduces sweating and bacteria formation**. The composition of the fabric allows easy sanitation and washing with quick drying times. The texture maintains the performance even at high loads and it has **self-supporting effect**. It is **fire retardant** according to reaction to fire Class 1 (non-flammable combustible materials) according to the MD of 26th June 1984 and subsequent amendments and additions.
 - » Removable cover is made by very breathable and high quality material both internally and externally. The inner padding has a **highly innovative leading edge structure**, that is three-dimensional, **lightweight**, providing extremely high air circulation, **optimal weight distribution** and excellent comfort. The external fabric is **highly resistant to abrasion** too. **Fire retardant** according to EN 1021-1:2014 and EN 1021-2:2014 standards.
- **Backrest** adjustable in height and tilt thanks to gas spring. Foldable for transport.
- Soft **armrests**, adjustable in height, tilt, width and rotation.
- **Multiaadjustable padded pelvic side supports** for greater containment. Adjustable in width, rotation and depth.
- **Safe adjustments without any risk of losing the screw.**
- **Seat** tilted by gas spring; adjustable in depth and width thanks to padded pelvic side supports.
- **Pelvic belt 947** provided as standard.
- **Legrest** adjustable in height and tilt with physiological fulcrum, such as to make the adjustment comfortable. Adjustable in depth with the seat.
- **Calf rests** adjustable in height and depth.
- **Footrests** tilting and adjustable in flexion-extension.
- **4 SWIVELLING WHEELS**, with the **NEW rubber coating** designed to prevent dirt from getting trapped in them and all equipped with parking brake. They allow you to move the Grillo Adaptive Seating effortlessly while the child is sitting.



Marked



Attention!
Read carefully the user and
maintenance manual

Please check off the interested products



MEASURES

Backrest upholstery colour:



Red



Pink




Blue

Frame colour: white

measures table (cm/in) with and without postural headrest




With the external headrest
863 don't raise the backrest
above the red line.



	A	B	C1	C2	D	E	F	G	L	M	Tot. weight	Max. Load	User height (H)
Medium	25-35 cm 9,8-13,7 in	28-43 cm 11-16,9 in	40-55 cm 15,7-21,6 in	40-70 cm 15,7-27,5 in	16-20 cm 6,2-7,8 in	20 cm 7,8 in	58 cm 22,8 in	74 cm 29,1 in	25-40 cm 9,8-15,7 in	70 cm 27,5 in	24,3 kg 53,5 lb	45 kg 99,2 lb	110-145 cm 43,3-57 in

↓ TO BE COMPLETED WITH THE HEADREST

☐ **865** Padded headrest adjustable in height.



☐ **835** Shaped headrest with lateral and occipital support adjustable in height.



☐ **863** Postural headrest That adapts to the conformation of the head. Adjustable in height, tilt, and depth.



TO FACILITATE THE DEVICE'S ASCENT AND DESCENT, CHOOSE ONE OF THE TWO GAS SPRINGS:

☐ **STANDARD (600N)** recommended for users up to 45 kg.

☐ **LIGHTER (400N)** recommended for users up to 35 kg.



If the headrest model and the gas spring are not selected, the high chair will be sent with the simpler 865 model and with the standard gas spring.



↓ ADDITIONAL COMPONENTS

947 Pelvic belt with variable angle (provided as standard)

☐ S ☐ M



920 Four point pelvic belt

☐ S ☐ M



838 Trunk side supports

Adjustable in height, width and transversally. Covered by breathable padding fabric, removable for washing.

☐ S ☐ M




☐ **868** Wrappable and flexible trunk supports

Adjustable in height, width and transversally. Covered by breathable padding fabric, removable for washing.



☐ **961** Elbow side supports.

Covered by breathable padding fabric, removable for washing.




853 Vest harness

☐ S



853 slim Four point shaped harness

☐ M




834R Adjustable abduction block in depth and abduction.

☐ mini/s ☐ M



834N Narrow abduction block

☐ mini/s ☐ M



☐ **824** Tray with wrap-around recess Adjustable in depth, height and tilt with armrests.



☐ **960** Heelrests
















☐ **962** Footrests + Heelrests



Notes: _____

C
O
M
P
O
N
E
N
T
S

M
E
A
S
U
R
E
S

<div>863</div> <div></div> <div>A=19 cm/7,4 in B=13-19 cm/5,1-7,4 in</div>	<div></div> <div>S/M A=18 cm/7 in</div>		
<div>947</div> <div></div> <div>small A=19 cm/7,4 in B=6 cm/2,3 in medium A=25,5 cm/10 in B=8 cm/3,1 in</div>	<div>920</div> <div></div> <div>small A=16 cm/6,2 in B=15,5 cm/6,1 in medium A=23 cm/9 in B=19 cm/7,4 in</div>	<div>838</div> <div></div> <div>small A=10 cm/3,9 in B=12 cm/4,7 in medium A=11 cm/4,3 in B=17 cm/6,6 in</div>	<div></div> <div>C=22-30 cm/8,6-11,8 in</div>
<div>868</div> <div></div> <div>small A=26 cm/10,2 in B=9 cm/3,5 in</div>	<div>853</div> <div></div> <div>small A=25 cm/9,8 in B=34 cm/13,3 in</div>	<div>853 slim</div> <div></div> <div>medium A=33 cm/12,9 in B=46 cm/18,1 in C=5 cm/1,9 in</div>	<div>834R</div> <div></div> <div>mini/s A=12 cm/4,7 in B=8-10 cm/3,1-3,9 in C=9 cm/3,5 in medium A=14 cm/5,5 in B=10-13 cm/3,9-5,1 in C=13 cm/5,1 in</div>
<div>834N</div> <div></div> <div>mini/s A=12 cm/4,7 in B=4 cm/1,5 in C=14 cm/5,5 in medium A=15 cm/5,9 in B=4 cm/1,5 in C=17 cm/6,6 in</div>	<div>824</div> <div></div> <div>A=16 cm/6,2 in B=38 cm/14,9 in C=53 cm/20,8 in D=64 cm/25,1 in</div>	<div>ARMRESTS ADJUSTABLE IN WIDTH AND ROTATION NEW</div> <div></div> <div>VIEW FROM ABOVE</div>	