

- BASE for OUTDOOR use
- BASE for INDOOR use

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER

MODULAR POSTURAL SYSTEM that can be positioned on OUTDOOR and INDOOR BASES



USER HEIGHT: 75-105 cm
(29,5-41,3 In)
MAX LOAD: 25 kg (55,1 lb)

FEATURES

Postural SEAT easily reversible with SERVO-ASSISTED TILTING AND RECLINING SYSTEM

- **Tilting seating unit** operated by gas spring, **removable** thanks to a practical, safe, and foolproof attachment/detachment system, quickly reversible to face the parent (outdoor base), **equipped with a visual safety system**. The backrest is adjustable in height, the seat in depth and the width can be adjusted via multi-adjustable lateral supports. All the adjustments are precise and continuous.
- The upholstery of the seat is extremely breathable, comfortable, durable, easy to wash, and quick to dry, consisting of:
 - » **Non-removable supporting upholstery:** the texture with an essential design allows an **excellent breathability**, reducing sweating and bacterial growth. The fabric composition allows an **easy sanitization**. **It maintains performance even at high loads and is fire retardant** according to class 1 fire reaction (non-flammable combustible materials) according to Ministerial Decree June 26, 1984, and subsequent amendments.
 - » **The removable cover** is made of breathable and **highly innovative** material both internally and externally. The internal fabric is three-dimensional, lightweight, and provides high air circulation. It allows an **optimal weight distribution and great comfort**. The external fabric is also highly resistant to abrasion. Fireproof according to EN 1021-1:2014 and EN 1021-2:2014 standards. Washable at 40°C. Easily removable for washing with **quick drying times**.
- **Reclining backrest** with gas spring and height-adjustable precisely and continuously. The backrest has a **restraining shape** and, thanks to the resistant bearing upholstery and **aluminum bars**, provides a comfortable and breathable support.
- **Integrated headrest** in the backrest to ensure the child always has support when the seating unit is tilted and/or the backrest is reclined.
- **Includes the mod. 947 pelvic belt with variable angle.**
- **Wide and padded armrests**, adjustable in height, width, and inclination (they follow the adjustments of the pelvis lateral supports).
- Adjustable **pelvic supports** for containment in width, rotation, and depth independently.
- Height and angle-adjustable **leg rest** with physiological pivot, making the adjustment comfortable. Adjustable in depth along with the seat.
- **Adjustable single footrest with heel stops** or **adjustable separate footrests** (component mod. 957).
- **Lateral shoulder supports** mod. 967 provided as standard in the outdoor version.



Attention!
Read carefully the user and maintenance manual

OUTDOOR BASE (CRASH TESTED according to ISO 7176/19 and ANSI RESNA WC/19 standards) **parent-facing reversible.**

- Compact and lightweight **frame** with a new innovative design, easily maneuverable and foldable, with a tilt-adjustable handlebar. **It reduces depth space.**
- Solid, lightweight, and smooth-running **wheels** on all surfaces: 17 cm (6,7 In) front wheels with direction locks, 25 cm (9,8 In) rear wheels, with shock absorbers and independent drum and parking brakes.

INDOOR BASE and LOWER INDOOR BASE adjustable in height.

- **Compact frame with height adjustment** - to allow approaching support and work surfaces - with gas spring adjustment while the child is seated. **The height adjustment system is internationally patented.**
- **Active tilt-in-space:** useful for positioning the child in an active position, to perform activities that stimulate straightening and activation of the posterior musculature.
- **Solid twin wheels** 7,5 cm (2,9 in), swiveling, **with rubber coating.** All equipped with parking brakes.

CHOOSE THE COLOR OF THE POSTURAL SEATING UNIT

☐ Coral pink



☐ Sky blue



☐ Brown



CHOOSE THE BASE (both can also be selected)

☐ Outdoor base




☐ Indoor base



☐ Lower indoor base




MEASURES in cm and (inches)




Seating unit on
outdoor base

A	B	C	D	E	F	G	H	I	L	M	N	O	P	Q	R	S	Frame weight	Seating unit weight	Total weight
34 (13,4)	16-26 (6,2-10,2)	57 (22,4)	42-52 (16,5-20,4)	51 (20)	72 (28,3)	45°	16-26 (6,2-10,2)	11,5-14 (4,5-5,5)	30°	15-25 (5,9-9,8)	31 (12,2)	65 (25,5)	69 (27,1)	49 (19,2)	57 (22,4)	80 (31,4)	10,8 kg (23,8 lb)	9,7 kg (21,3 lb)	20,5 kg (55 lb)



Seating unit on
indoor base

A	B	C	D	E	F	G	H	I	L	M	N	Frame weight	Seating unit weight	Total weight
34 (13,4)	16-26 (6,2-10,2)	53 (20,8)	42-52 (16,5-20,4)	45-65 (17,7-25,5)	65 (25,5)	45°	16-26 (6,2-10,2)	11,5-14 (4,5-5,5)	40°	15-25 (5,9-9,8)	40 (15,7)	14,3 kg (31,5 lb)	9,7 kg (21,3 lb)	24 kg (52,9 lb)



Seating unit on
lower indoor base

A	B	C *	D	E *	F	G	H	I	L	M	N	Frame weight	Seating unit weight	Total weight
34 (13,4)	16-26 (6,2-10,2)	60 (23,6)	42-52 (16,5-20,4)	36-48 (14-18,8)	63 (24,8)	45°	16-26 (6,2-10,2)	11,5-14 (4,5-5,5)	40°	15-25 (5,9-9,8)	26 (10,2)	9,8 kg (21,6 lb)	9,7 kg (21,3 lb)	19,5 kg (43 lb)

ADDITIONAL COMPONENTS for ALL BASES

☐ **863 Postural headrest**
adjustable in height, width, tilt and depth,
with adjustable side supports.
Size mini/S.



☐ **942 Ergonomic and mouldable headrest** adjustable in height, tilt and forward-backward with adjustable sides.
Size mini/S.




☐ **835 Padded shaped headrest** adjustable in height.



☐ **852 Headrest with parietal supports** adjustable in width.
Size mini/S.



☐ **838 Trunk side supports** adjustable in height, width and transversally.



967 Shoulders lateral supports. Standard on outdoor base.
☐ Optional component for the indoor base



☐ **868 Wrappable and flexible trunk supports** adjustable in height, width and transversally.



☐ **853 Vest harness**



☐ **853H Harness.**



947 pelvic belt with variable angle as standard.
Size mini.



☐ **834N Padded narrow abduction block.**
Size mini/S.



☐ **834R Adjustable abduction block.**
Size mini/S.



☐ **969 Postural pads kit.**



☐ **824 Transparent tray** adjustable in height, depth and tilt.



☐ **827 Foot straps** for single footrest.



☐ **957 Legrest with separate footrests.**
Adjustable in height, inclination and abduction independently.



☐ **960 Heelrests only with legrest 957.**



☐ **962 Footrests + Heelrests,** only with legrest 957.



☐ **892 Padded footplate covering,** for single footrest. Removable for washing.



☐ **972 Padded footrest height reducer** for single footrest.



ADDITIONAL COMPONENTS for OUTDOOR BASE

☐ **819 Versatile canopy.**



☐ **825 Rain cover.** Only with canopy 819.



☐ **839 Padded front handle** Adjustable in height and tilt.



☐ **818 Padded thermal cover.**



☐ **858 Shopping and vent tray basket.**



☐ **935 Oxygen bottle holder.**



☐ **891 Set of 4 tie down hooks** for the forward-facing transport of the pushchair in a motor vehicle (cars, buses...)



COMPLIANT WITH
ISO 7176/19 and
ANSI RESNA WC/19
STANDARDS



☐ **922** Bottle holder.



905 Hand brake.

☐ Lever Right

☐ Lever Left



ADDITIONAL COMPONENTS for INDOOR BASES

☐ **961** Elbow side supports.
Height adjustable.




☐ **824B**

Wide, resistant and hygienic table, suitable for food contact. Multiaadjustable.




863




mini/S A = 12 - 21 cm
4,7 - 8,2 in

942




mini/S A = 14 cm
5,5 in

852




mini/S A = 12 - 21 cm
4,7 - 8,2 in

835



mini A = 16 cm/6,3 in


838



mini
A = 9 cm/3,5 in
B = 7 cm/2,7 in


C = 16-24 cm
6,3 - 9,4 in

947




mini
A = 16 cm/6,3 in
B = 5 cm/1,9 in

868




mini
A = 23 cm/9 in
B = 9 cm/3,5 in

853H




A = 17-30 cm
6,7-11,8 in
B = 31 cm
12,2 in
C = 11,5-24 cm
4,5-9,4 in

853



mini
A = 20,5 cm/8 in
B = 28 cm/11 in


969



A = 19 cm/7,4 in
B = 7 cm/2,7 in
C = 2 cm/0,78 in


D = 7 cm/2,7 in
E = 30 cm/11,8 in
F = 2 cm/0,78 in

834R




mini/Small
A = 12 cm/4,7 in
B = 8-10 cm/3,1-3,9 in
C = 9 cm/3,5 in

834N




mini/Small
A = 12 cm/4,7 in
B = 4 cm/1,5 in
C = 14 cm/5,5 in

824B




A = 16 cm/6,2 in
B = 38 cm/14,9 in
C = 53 cm/20,8 in
D = 64 cm/25,1 in

824



A = 9 cm/3,5 in
B = 21 cm/8,2 in
C = 31 cm/12,2 in
D = 39 cm/15,3 in

858




A = 16 cm/6,2 in
B = 40 cm/15,7 in
C = 27 cm/10,6 in




5 kg./11 lb

957



A = 15° of abduction
on indoor base
A = 5° of abduction
on outdoor base

972



A = 11 - 21 cm
4,33-8,26 in