

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



Marked

Attention!
Read carefully the user and maintenance manual

FEATURES

- Compact and adjustable in height - to allow access to different table tops and work surfaces - thanks to the gas spring for the adjustment with the child seated. ONLY FOR INDOOR USE.
- **Frame** painted with non-toxic epoxy powders.
- **Safety system to avoid unintentional activation.**
- **Upholstery** extremely **breathable, resistant and easily washable**, even if not removable. **Fire retardant** according to EN 1021-1:2014 and EN 1021-2:2014.
- **Backrest** easily adjustable in height by sliding. Tilt with green levers.
- Soft **armrests** adjustable in height and tilt.
- **Seat** slightly tilted of 5° for better comfort. Adjustable in depth and width thanks to **958R multiadjustable padded pelvic side supports**.
- It can be **matched** to the corresponding measure of the **ORMESA Multifunction Table** (see its Order Form).

Please check off the interested products





MEASURES

Size: ☐ Small ☐ Medium

Upholstery colour: ☐ Pink ☐ Blue

Frame: Pebble grey



with 959 component,
set of 4 wheels with footrest

	A*	B	C	D	E	F	F1	G	G1	L	L1	User height (H)	Tot. weight	Max. Load
Small	30 cm 11,8 in	20-30 cm 7,8-11,8 in	30-45 cm 11,8-17,7 in	15-18 cm 5,9-7 in	36 cm 14,1 in	48 cm 18,8 in	50,5 cm 19,8 in	50 cm 19,6 in	52-64 cm 20,4-25,1 in	20-31 cm 7,8-12,2 in	28-39 cm 11-15,3 in	90-120 cm 35,4-47,2 in	10,7 kg 23,5 lb	35 kg 77,1 lb
Medium	35 cm 13,7 in	28-43 cm 11-16,9 in	35-55 cm 13,7-21,6 in	16-20 cm 6,2-7,8 in	41 cm 16,1 in	52 cm 20,4 in	54,5 cm 21,4 in	56,5 cm 22,2 in	62-73 cm 24,4-28,7 in	25-40 cm 9,8-15,7 in	33-48 cm 12,9-18,8 in	110-150 cm 43,3-59 in	13 kg 28,6 lb	45 kg 99,2 lb

* with component 958R multi-adjustable pelvic side supports: size SMALL 20-30 cm (7,8-11,8 in); size MEDIUM 25-35 cm (9,8-13,7 in)


ADDITIONAL COMPONENTS	947 Pelvic belt with variable angle ☐ mini ☐ S ☐ M	920 Four point pelvic belt ☐ S ☐ M	838 Trunk side supports Adjustable in height, width and rotation. ☐ mini ☐ S ☐ M	853 Vest harness ☐ mini ☐ S	853 slim Four point shaped harness ☐ M
	834N Narrow abduction block ☐ mini/s ☐ M	958R Padded pelvic side supports Adjustable in width, rotation and depth.	824 Tray with wrap-around recess Adjustable in depth, height and tilt with armrests.	823 Backrest extension for the head support	835 Shaped headrest with lateral and occipital support, adjustable in height. Only with 823 backrest extension.
	☐ 865 Padded headrest adjustable in height. Only with 823 backrest extension.	☐ 959 Set of 4 wheels with brakes as an alternative to the 4 black glides. 75 mm swivel wheels all equipped with brakes. Supplied with footrest.	⚠ Meet user's height and weight		

Notes: _____


COMPONENTS
MEASURES

947

mini A=16 cm B=5 cm A=6,2 in B=1,9 in	small A=19 cm B=6 cm A=7,4 in B=2,3 in	medium A=25,5 cm B=8 cm A=10 in B=3,1 in
--	---	---

920


small A=16 cm/6,2 in B=15,5 cm/6,1 in	medium A=23 cm/9 in B=19 cm/7,4 in
--	---


838

mini A=9 cm/3,5 in B=7 cm/2,7 in	small A=10 cm/3,9 in B=12 cm/4,7 in	medium A=11 cm/4,3 in B=17 cm/6,6 in
---	--	---


CHAIR size small
C= 18-26 cm/7-10,2 in

CHAIR size medium
C= 22-30 cm/8,6-11,8 in




853

mini A=20,5 cm/8 in B=28 cm/11 in	small A=25 cm/9,8 in B=34 cm/13,3 in
--	---

853 slim

medium A=33 cm/12,9 in B=46 cm/18,1 in C=5 cm/1,9 in
--


834N


mini/S A=12 cm/4,7 in B=4 cm/1,5 in C=14 cm/5,5 in	medium A=15 cm/5,9 in B=4 cm/1,5 in C=17 cm/6,6 in
--	--

823

small
A= 14 cm/5,5 in

medium
A= 20 cm/7,8 in



824

A= 16 cm/6,2 in	B= 38 cm/14,9 in	C= 53 cm/20,8 in	D= 64 cm/25,1 in
-----------------	------------------	------------------	------------------