

DISTRIBUTOR/AGENT	DATE OF ORDER
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



ISO 7176/19 and
ANSI RESNA WC/19



very breathable
and fire retardant
upholstery



Marked



Attention!
Read carefully the user and
maintenance manual

FEATURES

- **More sizes in one frame**
- **Frame:** light and umbrella folding system, in aluminium. It is painted with non-toxic epoxy powders. When folded Trollino is compact and it allows you to carry it as if it was a travel trolley.
- **Push handles** with **manual and parking brakes**.
- Resistant and durable **backrest** guaranteed by the cross on the back. Reclinable on three positions and easily adjustable in height.
- Slightly tilted **seat**. It is easily adjustable in depth and width thanks to padded pelvic side supports (additional components 958). The cross below ensures stability, comfort and durability.
- Extremely breathable **upholstery** (not removable). It is resistant and easily washable. **Fire Retardant** according to *EN 1021-1:2014* and *EN 1021-2:2014*.
- User manual **bag** and wrench on the back of the backrest.
- **Legrest** adjustable in height.
- Elevatable and swivel **footrests**, without encumbrance.
- **Solid wheels:** 17 cm (7 in) front swiveling wheels and 25 cm (10 in) rear wheels with **independent drum brakes**. **Overcoming step bar** on both rear wheels.

Please check off the interested products



MEASURES


Backrest colour: Seat colour: **Black** Frame colour: **Black**

- ☐ Red
☐ Pink
☐ Blue
☐ Black

Unit of Measurement: cm (in)



Total Weight: 14,8 kg (32,6 lb)
Max. Load: 55 kg (121 lb)



A: 36 cm (14 in). For different width (32 and 28 cm - 12,6 and 11 in) please select 958 padded pelvic side supports.
*Height adjustment range with add-on 951 max. 105 cm (41,3 in).

☐ **951** Height adjustable push handles. Tools free!
COMPONENT TO BE ORDERED WITH THE PRODUCT



COMPONENTS TO CHOOSE ALSO AFTERWARDS

835 Shaped and
enveloping headrest
adjustable in height.

☐ S-M




865 Padded
headrest adjustable in
height.

☐ S-M




838 Trunk side
supports adjustable in
height, width and
inclination.

☐ mini ☐ S



868 Wrappable and
flexible trunk support,
multiadjustable in height,
width and inclination.

☐ S



☐ **958** Padded
pelvic side supports
32 cm (12,6 in) seat width.

☐ **Extra Padding**
28 cm (11 in) seat width




947 Pelvic belt with
variable angle.

☐ S ☐ M




920 4 point pelvic.
belt

☐ S ☐ M




853 Vest harness.

☐ mini ☐ S




853 slim
4 point shaped
harness.


☐ M



☐ **834N** Padded
and removable
abduction block.



☐ **839**
Front handlebar
adjustable in tilt.



☐ **960**
Heel rests.



☐ **956** Footrests
adjustable in tilt.



☐ **858**
Shopping basket.



☐ **819** Elevating
Canopy with bag.
Foldable with the
frame.



☐ **825** Rain cover
to fix on the 819 canopy.



☐ **825B** Back
and lateral protection
clothe to fix on the 819
canopy.




☐ **818** Thermic
cover.



☐ **968** Padded
covers for armrests.
Easily removable for
washing.



☐ **844** Padded
side protection
covers.



☐ **812** Front
wheel direction
locks.



☐ **891** Set of 4 tie down hooks
for the forward-facing transport of the pushchair in a
motor vehicle.



COMPONENTS MEASURES

Notes.....