

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

FEATURES

- **FRAME** made of chromed steel, painted with non-toxic epoxy powder paint. Narrow base to facilitate entering elevators.
- **The seat and backrest are covered with visco-elastic fire-proof material. It is resistant, breathable and anti-slip. It is made of a special and exclusive ORMESA formula:**
  - 1) It adapts to the user's body shapes and it prevents its sliding, with the result of an **EXCELLENT COMFORT and KEEPING POSTURE**, even when changing the backrest recline and tilt (it adapts to the user's body shapes maintaining its support function).
  - 2) Easily washable to keep the **HYGIENE** of the user; it lets **AIR PASS** keeping the **OPTIMAL MICROCLIMATE** that **HELPS PREVENTING BED SORES**.
  - 3) Durable and fireproof, **GUARANTEE OF DURABILITY and SAFETY**.
- Ergonomic, breathable, non-toxic, fireproof **UPHOLSTERY**; It does not interfere with the benefits of the visco-elastic cover. The upholstery is fixed with buttons, easy to remove for machine washing.
- Ergonomic **BACKREST**, it holds the shapes of the user. **The special structure of the backrest gives a proper alignment to the trunk and it doesn't generate any pressure on the spine.** Adjustable in tilt by levers with safety pin. The pin allows blocking the levers to maintain the proper adjustment for the user, height adjustable, folding for transportation.
- **Low SEAT** and slightly inclined. Adjustable in depth on four positions.
- Wide and wrap **HEADREST**, padded, multi adjustable (in height, tilt, depth, lateral traslation); thus assuring the **proper alignment between head and trunk**, without interfering with the user's view. Anti-bacterial and Hygienic.
- **ARMRESTS**, padded, easily and stably adjustable in height thanks to the lever.
- **PUSH HANDLE** adjustable in tilt/ height, it helps the care.
- **CONTINUOUS AND WITHOUT EFFORT ADJUSTMENTS:** tilting, backrest recline and legrest tilt are operated by independent gas springs.
- **SIDES** with scissor opening. Thanks to the special joint, the armrests are positioned behind the backrest, eliminating completely all obstacles to the transfer.
- **LEG RESTS:** ergonomic, removable, turnable and adjustable in tilt without effort, activated with gas springs.
- **CALFRESTS** adjustable in depth and free tilt to accommodate the movements of the leg.
- **FOOT RESTS**, adjustable in flexion-extension and height.
- **8" front solid swiveling wheels and 12" rear wheels;** drum brakes are operated by a single control pedal.



Marked

Tie down hooks (4 rings) for the front facing transport of the user in moving vehicles (private cars, conformed buses etc.) **AS STANDARD**



Please check off the interested products:




**Attention!**  
Read carefully the user and maintenance manual

SIZE  36  40  45  50

UPHOLSTERY COLOUR:  dark grey  light grey

FRAME COLOUR:  white  black

MEASURES







\*footrest adjusted to the maximum length

	A seat bottom width	B seat bottom depth	C backrest width	D overall width	H seat bottom inclination-tilt	I backrest inclination	Total weight	Max. Capacity	User height
<b>size 36</b>	36 cm 14,1 in	40-43-46-49 cm 15,7-16,9-18,1-19,2 in	47 cm 18,5 in	62 cm 24,4 in	from 10° to 35°	from 10° to 38°	38 kg 83,7 lb	120 kg 264,5 lb	from 140 to 190 cm from 55,1 to 74,8 in
<b>size 40</b>	40 cm 15,7 in	40-43-46-49 cm 15,7-16,9-18,1-19,2 in	47 cm 18,5 in	62 cm 24,4 in	from 10° to 35°	from 10° to 38°	38 kg 83,7 lb	120 kg 264,5 lb	from 140 to 190 cm from 55,1 to 74,8 in
<b>size 45</b>	45 cm 17,7 in	40-43-46-49 cm 15,7-16,9-18,1-19,2 in	52 cm 20,4 in	67 cm 26,3 in	from 10° to 35°	from 10° to 38°	41 kg 90,3 lb	130 kg 286,6 lb	from 140 to 190 cm from 55,1 to 74,8 in
<b>size 50</b>	50 cm 19,6 in	40-43-46-49 cm 15,7-16,9-18,1-19,2 in	57 cm 22,4 in	72 cm 28,3 in	from 10° to 35°	from 10° to 38°	42 kg 92,5 lb	140 kg 308,6 lb	from 140 to 190 cm from 55,1 to 74,8 in

Hyphen separated numbers indicate a minimum and maximum range of adjustment

**ADDITIONAL COMPONENTS**

<p><input type="checkbox"/> <b>824-R</b>                  Reclinable tray, adjustable in tilt, backwards - forwards, with security block</p> <p>It is made by resistant and hygienic material, design with no edges, it doesn't keep residual food.</p>  <p><i>igienico - hygienic                  resistente e leggero                  durable and lightweight                  ergonomico - ergonomic                  regolazioni semplici e sicure                  simple and safe adjustments</i>  <b>New design!</b></p>	<p><input type="checkbox"/> <b>894</b> 45° Pelvic belt</p> 	<p><input type="checkbox"/> <b>834</b> Padded abduction block</p> 	<p><input type="checkbox"/> <b>892</b> Footplate covering</p> 	
<p><input type="checkbox"/> <b>905</b> Set of hand brakes with single lever</p> 	<p><input type="checkbox"/> <b>914</b> Footrest lock</p> 	<p><input type="checkbox"/> <b>896</b> Tray covering</p> 	<p><input type="checkbox"/> <b>916</b> Padded leg panel</p> 	<p><input type="checkbox"/> <b>932</b> Anti-tip device</p> 
<p><input type="checkbox"/> <b>933</b> IV Pole</p> 	<p><input type="checkbox"/> <b>936</b> Headrest cover</p> 	<p><input type="checkbox"/> <b>938</b> Water-proof cover</p> 	<p><b>975</b> Multi adjustable side-temporo support. It can be positioned on the right or left side.  <input type="checkbox"/> 1 pc. <input type="checkbox"/> 2 pc.</p> <p><b>NEW</b></p> 	<p><input type="checkbox"/> <b>976</b> Multi adjustable side-temporo-fronto support. It can be positioned on the right or left side.</p> <p><b>NEW</b></p> 

Notes: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_