

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



according to the
ISO 7176/19

NEW



CE marked



Attention!
Read carefully the user and
maintenance manual

FUNCTIONS


- Lightweight aluminium **frame**, painted with non-toxic epoxy powders. It **folds just like an umbrella**.
- **Semi-rigid backrest**, reclines in three positions with zippers.
- **Height adjustable Handles, Tools free.**
- **Semi-rigid seat** bottom, slightly inclined.
- **Padded upholstery** removable for washing with **removable aluminium reinforcement ribs**. Extremely breathable and quick to dry, excellent for **load distribution**; **Fireproof fabric** according to EN 1021-1:2014 and EN 1021-2:2014.
- **Net basket** on the back of the backrest.
- Removable **footrest** height adjustable.
- **Legrest strap**.
- Puncture-proof **wheels**: front 15 cm swivelling, rear 20 cm with independent brakes **and overcoming step bar**.

Please check off the interested products:

Upholstery colour: **Black and red** Frame colour: **Black**


- ☐ **CLIP 1**
- ☐ **CLIP 2**
- ☐ **CLIP 3**

Measures in cm (Inches):



Seat, top view

17 cm / 6,7 in





32 cm / 12,6 in

Seat width with padded sides,
removable with zip.



NEW



SIZE	Max. Load	Total Weight	A	B	B1*	C	D	E	F	G	H	I	L	M	N	O	P	Q
1	30 kg 66,1 lb	10,3 kg 22,7 lb	40 (15,7)	51 (20)	39 (15,3)	101-108 (39,7-42,5)	0°-12°	25 (9,8)	55 (21,6)	18°	46 (18,1)	12-20 (4,7-7,8)	78 (30,7)	84 (33)	114 (44,8)	30 (11,8)	34 (13,3)	110-125 (43,3-49,2)
2	35 kg 77,1 lb	10,5 kg 23,1 lb	40 (15,7)	51 (20)	39 (15,3)	101-108 (39,7-42,5)	0°-12°	30 (11,8)	59 (23,2)	18°	48 (18,1)	21-31 (8,2-12,2)	78 (30,7)	92 (36,2)	114 (44,8)	30 (11,8)	34 (13,3)	120-140 (47,2-55,1)
3	40 kg 88,1 lb	10,7 kg 23,5 lb	40 (15,7)	51 (20)	39 (15,3)	106-113 (41,7-44,4)	0°-12°	35 (13,7)	63 (24,8)	18°	49 (19,2)	25-35 (9,8-13,7)	84 (33)	98 (38,5)	118 (46,4)	34 (13,3)	36 (14)	135-155 (53,1-61)

* **Seat width**, without padded sides (B1)

☐ **835** Headrest with parietal supports.



☐ **816** Trunk lateral containment pads.




☐ **853** Vest harness.
(S for CLIP 1 and 2 - M for CLIP 3)



☐ **906** Five point harness.
(S for CLIP 1 and 2 - M for CLIP 3)



☐ **903** Five point vest harness.
(S for CLIP 1 and 2 - M for CLIP 3)



☐ **894** 45° Pelvic belt.
(S for CLIP 1 and 2 - M for CLIP 3)



☐ **834N** Padded and removable abduction block.
☐ mini/S ☐ M



☐ **858** Shopping basket.



☐ **827** Foot straps.



☐ **812** Front wheel direction locks.



☐ **819** Canopy.



☐ **818** Thermal cover.




☐ **825** Rain cover.
(to fix to the canopy 819)



☐ **891** Set of 4 tie down hooks for the forward-facing transport of the pushchair in a motor vehicle (private vans/cars, buses..)




853



Small
A=30 cm/11,8 in
B=31,5 cm/12,4 in

Medium
A=39,5 cm/15,5 in
B=39 cm/15,3 in

834N



mini/S
A=12 cm/4,7 in
B=4 cm/1,5 in
C=14 cm/5,5 in

Medium
A=15 cm/5,9 in
B=4 cm/1,5 in
C=17 cm/6,6 in

Notes.....

.....

.....

.....

.....