

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



Marked

FEATURES

- Compact and adjustable in height - to allow access to different table tops and work surfaces - thanks to the gas spring for the adjustment with the child seated. ONLY FOR INDOOR USE.
- **Frame** painted with non-toxic epoxy powders.
- **Safety system to avoid unintentional activation.**
- **Upholstery** extremely **breathable, resistant and easily washable**, even if not removable. **Fire retardant** according to EN 1021-1:2014 and EN 1021-2:2014.
- **Backrest** easily adjustable in height by sliding. Tilt with green levers.
- Soft **armrests** adjustable in height and tilt.
- **Seat** slightly tilted of 5° for better comfort. Adjustable in depth and width thanks to **958R multiadjustable padded pelvic side supports**.
- It can be **matched** to the corresponding measure of the **ORMESA Multifunction Table** (see its Order Form).



Attention!
Read carefully the user and maintenance manual

Please check off the interested products





MEASURES


Size: ☐ Small ☐ Medium

Upholstery colour: ☐ Pink ☐ Blue

Frame: Pebble grey



with 959 component,
set of 4 wheels with footrest




	A*	B	C	D	E	F	F1	G	G1	L	L1	User height (H)	Tot. weight	Max. Load
Small	30 cm 11,8 in	20-30 cm 7,8-11,8 in	30-45 cm 11,8-17,7 in	15-18 cm 5,9-7 in	36 cm 14,1 in	48 cm 18,8 in	50,5 cm 19,8 in	50 cm 19,6 in	52-64 cm 20,4-25,1 in	20-31 cm 7,8-12,2 in	28-39 cm 11-15,3 in	90-120 cm 35,4-47,2 in	10,7 kg 23,5 lb	35 kg 77,1 lb
Medium	35 cm 13,7 in	28-43 cm 11-16,9 in	35-55 cm 13,7-21,6 in	16-20 cm 6,2-7,8 in	41 cm 16,1 in	52 cm 20,4 in	54,5 cm 21,4 in	56,5 cm 22,2 in	62-73 cm 24,4-28,7 in	25-40 cm 9,8-15,7 in	33-48 cm 12,9-18,8 in	110-145 cm 43,3-57 in	13 kg 28,6 lb	45 kg 99,2 lb

* with component 958R multi-adjustable pelvic side supports: size SMALL 20-30 cm (7,8-11,8 in); size MEDIUM 25-35 cm (9,8-13,7 in)

ADDITIONAL COMPONENTS	947 Pelvic belt with variable angle ☐ mini ☐ S ☐ M	920 Four point pelvic belt ☐ S ☐ M	838 Trunk side supports Adjustable in height, width and rotation. ☐ mini ☐ S ☐ M	853 Vest harness ☐ mini ☐ S	853 slim Four point shaped harness ☐ M
	834N Narrow abduction block ☐ mini/s ☐ M	☐ 958R Padded pelvic side supports Adjustable in width, rotation and depth.	☐ 824 Tray with wrap-around recess Adjustable in depth, height and tilt with armrests.	☐ 823 Backrest extension for the head support	☐ 835 Shaped headrest with lateral and occipital support, adjustable in height. Only with 823 backrest extension.
	☐ 865 Padded headrest adjustable in height. Only with 823 backrest extension.	☐ 959 Set of 4 wheels with brakes as an alternative to the 4 black glides. 75 mm swivel wheels all equipped with brakes. Supplied with footrest.	⚠ Meet user's height and weight		


Notes: _____

947




mini	small	medium
A=16 cm	A=19 cm	A=25,5 cm
B=5 cm	B=6 cm	B=8 cm
A=6,2 in	A=7,4 in	A=10 in
B=1,9 in	B=2,3 in	B=3,1 in

920




small
A=16 cm/6,2 in
B=15,5 cm/6,1 in
medium
A=23 cm/9 in
B=19 cm/7,4 in

838




mini
A=9 cm/3,5 in
B=7 cm/2,7 in
small
A=10 cm/3,9 in
B=12 cm/4,7 in
medium
A=11 cm/4,3 in
B=17 cm/6,6 in

CHAIR size small
C= 18-26 cm/7-10,2 in
CHAIR size medium
C= 22-30 cm/8,6-11,8 in




853




mini
A=20,5 cm/8 in
B=28 cm/11 in
small
A=25 cm/9,8 in
B=34 cm/13,3 in

853 slim



medium
A=33 cm/12,9 in
B=46 cm/18,1 in
C=5 cm/1,9 in


834N




mini/S
A=12 cm/4,7 in
B=4 cm/1,5 in
C=14 cm/5,5 in
medium
A=15 cm/5,9 in
B=4 cm/1,5 in
C=17 cm/6,6 in

823

small
A= 14 cm/5,5 in
medium
A= 20 cm/7,8 in



824



A= 16 cm/6,2 in
B= 38 cm/14,9 in
C= 53 cm/20,8 in
D= 64 cm/25,1 in