



• MADE IN ITALY SINCE 1980 •



ORDER FORM

Frame runner
Assistive Device for running

Size L

ALUMINIUM

NEW

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:

FEATURES

- **Aluminum lightweight frame and fork.**
- **Adjustable handlebar** in height, backward-forward and reversible.
- **Ergonomic handles**, to grant grip and comfort.
- **Steering limiters with springs**, that allow the handlebar to return to the central position. It can be disengaged to allow more steering. It's possible to choose between two different elastic strengths.
- **Memory foam padded trunk support**, adjustable in height, backward-forward and free-angle tilt, with breathable padded upholstery, easily removable for washing.
- **Front disc brake with locking device** for ascent and descent, that can be operated by the lever either on the right or on the left side of the handlebar. The distance between the brake lever and the handlebar is adjustable by screw.
- **Sports saddle** adjustable in height, tilt and anterior-posterior direction.
- **Pneumatic road-type wheelset 28"** with reflector stripes, easily removable with key.



Product marked **CE**

WHAT TO WEAR





Attention!
Read carefully the user and maintenance manual

Please check off the interested products: ☒

☐ Size L

measures in cm and (inches)

MEASURES



Size	A	B	C	D	E	F	G	H	Frame weight	Total weight without wheels	Max Load
L	62 (24,4)	92 (36,2)	95-107 (37,4-42,1)	192 (75,6)	71-84 (27,9-33)	72,5-92 (28,5-36,2)	127 (50)	150-190 (59-74,8)	5,7 kg (12,5 lbs)	10,2 kg (22,5 lbs)	100 kg (220 lbs)



OPTIONAL COMPONENTS, ONLY CHOOSABLE WHEN ORDERING THE FRAME RUNNER

<p>Steering limiter spring</p> <p>☐ Soft</p>  <p>☐ Hard</p>  <p>If no specification is provided, the Frame Runner will be supplied with the hard spring fitted as standard.</p>	<p>Please choose the multiadjustable, padded trunk support size.</p> <p>☐ M ☐ L</p>  <p>If no size is specified, size L will be provided.</p>	<p>Please choose the handlebar size.</p> <p>☐ M ☐ L</p>  <p>Please choose the position of the brake on the handlebar or on the ergonomic handles (880).</p> <p>☐ Left ☐ Right</p> <p>In absence of choice, the brake will be mounted on the right.</p>	<p>Please choose the seatpost length.</p> <p>** Does not cover the entire adjustment range.</p> <p>☐ **15 cm (5,9")</p> <p>☐ **24 cm (9,5")</p> <p>☐ 28 cm (11")</p>  <p>In absence of choice, the 28 cm (11") seatpost will be supplied.</p>
--	--	--	--

<p>☐ 880 Ergonomic handles adjustable in inclination independently (1 pair). Same size as the Handlebar.</p> <p>Choose the position of the brake on Ergonomic Handles.</p> 	<p>875 Wrist retainer with support and multi-adjustable glove independently.</p> <p>☐ Left ☐ Right</p> 	<p>845 Forearm support.</p> <p>☐ Left ☐ Right</p> 	<p>102 Trunk support wide strap, width adjustable; to be applied to the trunk support.</p> <p>☐ S-M ☐ L</p> 
<p>☐ 876 Monocycle saddle.</p> 	<p>☐ 103 Monocycle saddle padding.</p> 	<p>☐ 114 Comfort saddle.</p> 	<p>☐ 113 Comfort saddle cover.</p> 
<p>☐ 104 Frame side paddings.</p> 	<p>☐ 105 Frame rear paddings.</p> 	<p>☐ 106 Under saddle padding.</p> 	<p>☐ 107 Parking brakes, right and left.</p> 

☐ **108** Front led lamp.
(with USB charging)



☐ **109** Back led lamp.
(with USB charging)




☐ **110** Bottle holder.
Water bottle not included.




☐ **111** Transparent rear wheels
spoke protectors, right and left.
Applicable to standard wheels
only.



☐ **117** Transport bag with handles, integrated wheels and full side opening with zipper, for Frame Runner with wheels removed.



A	B	C
74 cm/29 in	132 cm/51,9 in	75 cm/29,5 in




It can be folded for storage.

☐ **118** Set of 2 bags for front and
rear wheels, with zipper.




A	B
80/31,5 in	87 cm/34,5 in

Padded trunk support




	A	B	C
M	16 cm/6,3 in	32 cm/12,6 in	5,8 cm/2,3 in
L	20 cm/7,9 in	32 cm/12,6 in	5,8 cm/2,3 in

handlebar




	A
M	55 cm/21,6 in
L	64 cm/25,2 in

102




	A	B
S-M	8 cm/3,14 in	45-59 cm/17,7-23,2 in
L	8 cm/3,14 in	63,5-84 cm/25-33 in

876



A	B	C
13,5 cm/5,3 in	31 cm/12,2 in	9 cm/3,5 in

114



A	B
25,5 cm/10 in	20 cm/7,8 in