

DISTRIBUTOR/AGENT:	DATE OF ORDER:
	DATE OF DELIVERY:
	NAME OF THE ADDRESSEE/USER:



Product marked



Attention!
Read carefully the user and
maintenance manual

FEATURES

- **Adjustable handlebar** in height, backward-forward and reversible.
- **Wide and ergonomic handles**, to increase comfort.
- **Steering limiters with springs**, that allow the handlebar to return to the central position. Releasing to allow more steering. It's possible to choose between two different elastic strengths.
- **Memory foam padded trunk support**, adjustable in height, backward-forward and free-angle tilt, with breathable padded upholstery, easily removable for washing.
- **Front brake with locking device** for ascent and descent, that can be operated by the lever either on the right or on the left side of the handlebar. The distance between the brake lever and the handlebar is adjustable by screw.
- **Ergonomic saddle** adjustable in height, tilt and anterior-posterior direction.
- **Pneumatic road-type or MTB wheelset** with reflector stripes.
Quick release front wheel.
Rear wheels easily removable with key.


Please check off the interested products: ☒ _____

MEASURES

☐ size M

☐ size L

measures in cm and (inches)



	A	B	C	D ¹	E	F ²	G	H	Weight ³	Max Load
M	55 (21,6)	86 (33,8)	84-97 (33-38,1)	159 (62,6)	53-64 (20,8-25,1)	58-73 (22,8-28,7)	106 (41,7)	120-150 (47,2-59)	17 kg (37,4 lbs)	70 kg (154 lbs)
L	62 (24,4)	92 (36,2)	95-107 (37,4-42,1)	187 (73,6)	71-84 (27,9-33)	70-92 (27,5-36,2)	127 (50)	150-190 (59-74,8)	18 kg (39,7 lbs)	100 kg (220 lbs)

Dimensions measured with road-type wheelset.

¹ Size M with 26" wheel, D_{max}=164 cm (64,5") - Size L with 28" wheel, D_{max}=192 cm (75,5").

² With standard saddle. It is possible to further adjust the height position of the saddle (±3 cm - 1,18") by acting on the fixing clamp.

³ Weight calculated without saddle abatement system.

OPTIONAL COMPONENTS, ONLY CHOOSABLE WHEN ORDERING THE FRAME RUNNER

☐ Saddle abatement system , to facilitate the ascent and descent from the aid. ☐ Control lever oriented downward		  <p>NEW</p>	<p>Please choose the multiadjustable, padded trunk support size.</p> <p>☐ S ☐ M ☐ L</p>  <p>If no size is selected, the corresponding size will be supplied.</p>
<p>Standard configuration: control lever oriented upward (as shown in the picture). Select the option above to request downward lever orientation.</p>			
<p>Please choose the handlebar size.</p> <p>☐ S ☐ M ☐ L</p>  <p>Please choose the position of the brake on the handlebar or on the ergonomic handles (880).</p> <p>☐ Left ☐ Right</p> <p>In absence of choice, the brake will be mounted on the right.</p>	<p>Please choose the seatpost length.</p> <p>** The 15 cm (5,9") seatpost does not cover the entire adjustment range.</p> <p>☐ **15 cm (5,9") ☐ 28 cm (11")</p>  <p>In absence of choice, the 28 cm (11") seatpost will be supplied.</p>	<p>Steering limiter spring</p> <p>☐ Soft</p>  <p>☐ Hard</p>  <p>NEW</p> <p>In absence of choice, size M will be supplied with the soft spring and size L with the hard spring.</p>	

OPTIONAL COMPONENTS, ALSO CHOOSABLE SEPARATELY

101/101MTB Wheels dimensions: POSSIBLE CONFIGURATIONS

☐ **MTB tires** (In absence of choice, Frame Runner will be supplied with road-type wheelset).

Frame runner M

In absence of choice, the M size will be supplied with the trio of 24" road-type wheels.


- ☐ Front wheel with fork 20" - Rear wheels 24"
- ☐ Front wheel with fork 24" - Rear wheels 26"

Frame runner L

In absence of choice, the L size will be supplied with the trio of 26" road-type wheels.


- ☐ Front wheel with fork 24" - Rear wheels 24"
- ☐ Front wheel with fork 24" - Rear wheels 26"
- ☐ Front wheel with fork 26" - Rear wheels 24"
- ☐ Front wheel with fork 26" - Rear wheels 28"
- ☐ Front wheel with fork 28" - Rear wheels 28"

OPTIONAL COMPONENTS, ALSO CHOOSABLE SEPARATELY


☐ 880 Ergonomic handles adjustable in inclination independently (1 pair). Same size as the Handlebar. Choose the position of the brake on Ergonomic Handles. 	☐ 875 Wrist retainer with support and multi-adjustable glove independently. <div> ☐ Left ☐ Right </div> 	☐ 845 Forearm support. <div> ☐ Left ☐ Right </div> 	☐ 102 Trunk support wide strap, width adjustable; to be applied to the trunk support. <div> ☐ S-M ☐ L </div> 
☐ 876 Monocycle saddle. 	☐ 103 Monocycle saddle padding. 	☐ 113 Standard saddle cover. 	☐ 115 Forward/backward adjustable seatpost. Allows the saddle position to be adjusted forward or backward according to the user's needs. <div> ☐ M Only for Frame Runner M ☐ L Only for Frame Runner L NEW </div> 
☐ 104 Frame side paddings. 	☐ 105 Frame rear paddings. 	☐ 106 Under saddle padding. 	☐ 107 Parking brakes, right and left. ☐ 107MTB Parking brakes, right and left. (MTB WHEELS ONLY) 
☐ 108 Front led lamp. (with USB charging) 	☐ 109 Back led lamp. (with USB charging) 	☐ 110 Bottle holder. Water bottle not included. 	☐ 111 Transparent rear wheels spoke protectors, right and left. 

OPTIONAL COMPONENTS, ALSO CHOOSABLE SEPARATELY

- ☐ **117** Transport bag with handles, integrated wheels and full side opening with zipper, for Frame Runner with wheels removed.



A	B	C
74 cm/29 in	132 cm/51,9 in	75 cm/29,5 in



It can be folded for storage.

- ☐ **118** Set of 2 bags for front and rear wheels, with zipper.




A	B
80/31,5 in	87 cm/34,5 in

Notes:-----


COMPONENT MEASUREMENTS

Padded trunk support




	A	B	C
S	12 cm/4,7 in	24 cm/9,4 in	3,6 cm/1,4 in
M	16 cm/6,3 in	32 cm/12,6 in	5,8 cm/2,3 in
L	20 cm/7,9 in	32 cm/12,6 in	5,8 cm/2,3 in

Handlebar



	A
S	50 cm/19,6 in
M	55 cm/21,6 in
L	64 cm/25,2 in


115



	A	B
M	20 cm/7,87 in	20,5 cm/8,07 in
L	13 cm/5,12 in	20,5 cm/8,07 in


Saddle depth adjustment range:
M: 6 cm/2,36 in
L: 10 cm/3,94 in

102



	A	B
S-M	8 cm/3,14 in	45-59 cm/17,7-23,2 in
L	8 cm/3,14 in	63,5-84 cm/25-33 in

876



A	B	C
13,5 cm/5,3 in	31 cm/12,2 in	9 cm/3,5 in